

BD SmartSite™ Extension Set

SKU/REF 20019E7D
SKU/REF Name BD SmartSite™ 2-Way Extension Set - 24 cm



Specifications

Total Length (cm)	24
Priming Volume (mL)	0.6
Tube Inner Diameter (mm)	1
Ways	2-Way
Clamp	2
Luer Lock	Rotating Male
Needle Free Connector	BD SmartSite™
Pressure Rating	325 PSI @ 10 mL per second
PVC	Yes
DEHP	DEHP-Free
Shelf Life	3 Years
Sterilization Method	Gamma Irradiation
Country Of Manufacture	Mexico
Classification	Ila
Quantity per Box	100

CE 2797

<https://www.bd.com/en-eu/bd-smartsite-extension-set/20019e7d?return=true>

Becton Dickinson Sverige AB, Marieviksgatan 25, 11743 Stockholm
bdsweden@bd.com

bd.com

BD, BD-logo, NeutroX, NeutraClear är varumärken som tillhör Becton, Dickinson and Company eller dess dotterbolag
© 2022 BD. Alla rättigheter förbehållna. BD-59783 March 2022

