


IV3000 product range

Selection of the correct dressing for insertion sites is advised as a best practice to prevent infections¹. IV3000 offers an extensive range of dressings designed to meet the most common IV therapy practices. Our range includes four types of dressings in various sizes:

Ported dressings

A transparent film dressing with moisture-responsive technology^{2,3} that incorporates non-woven wings which reinforce fixation around the catheter site. The dressing has an opening designed to host the cannula's port. The dressing also includes two sterile securing strips for catheter or tube fixation and a sterile documentation label.^{4,5}




Sizes*	S&N code	Quantity per carton	Recommended indication
5cm x 6cm	66004011	100	Paediatric
7cm x 9cm	4006	100	Peripheral
9cm x 12cm (oval)	66004009	50	Central
11cm x 14cm	66800512	25	PICC

*Availability and product code could vary per country

Non-ported dressings

A transparent film dressing with moisture-responsive technology^{2,3} with two securing strips and documentation label.^{4,5}




Sizes*	S&N code	Quantity per carton	Recommended indication
5cm x 4cm**	4654	100	Paediatric
6cm x 7cm	4007	100	Peripheral
10cm x 12cm	4008	50	Central

*Availability and product code could vary per country

**This product size does not include securing strips and documentation label

Frame Delivery

A transparent film dressing with moisture-responsive technology^{2,3}. The frame delivery system allows user to position the dressing over the catheter with immediate visualization of the IV site. The dressing includes two securing strips and documentation label.




Sizes*	S&N code	Quantity per carton	Recommended indication
6cm x 7cm	59410082	100	Peripheral
10cm x 12cm	59410882	50	Central

*Availability and product code could vary per country

Reinforced handle (orange-handles)

A transparent film dressing with moisture-responsive technology^{2,3} with edges of orange tape. These can either be torn from the edges of the dressing and used for extra security of the catheter or left in place as part of the dressing.



Sizes*	S&N code	Quantity per carton	Recommended indication
6cm x 8cm	4923	100	Peripheral
6cm x 8.5cm	4924	100	Peripheral
10cm x 14cm	4925	10	Central
10cm x 14cm	4973	50	Central
10cm x 20cm	4649	50	PICC/Epidural

*Availability and product code could vary per country

References