


Codman® Surgical Patties and Strips

*The benchmark
in quality*



CODMAN Surgical Patties, sterile, X-ray detectable, 10 patties per package, 20 packages per dispenser box.

80-1396	5/16"	(7mm)
80-1397	1/4" x 1 1/2"	(7mm x 38mm)
80-1398	1/4" x 3"	(7mm x 76mm)
80-1399	1/4" x 1/4"	(7mm x 7mm)
80-1400	1/2" x 1/2"	(13mm x 13mm)
80-1401	3/4" x 3/4"	(19mm x 19mm)
80-1402	1/2" x 1"	(13mm x 25mm)
80-1403	1" x 1"	(25mm x 25mm)
80-1404	1/2" x 1 1/2"	(13mm x 38mm)
80-1406	1/2" x 2"	(13mm x 52mm)
80-1407	1/2" x 3"	(13mm x 76mm)
80-1408	1" x 3"	(25mm x 76mm)
80-1409	3" x 3"	(76mm x 76mm)

CODMAN® Surgical Patties and Strips

CODMAN Surgical Patties and Strips provide the neurosurgeon with softness, gentleness, protection and absorbency to meet the demands of delicate intracranial and spinal surgery.

Absorbency

CODMAN Patties and Strips absorb more than 5 times their weight in less than a second.

Precision Cut

CODMAN Patties and Strips' precision cut provides smooth edges and helps prevent linting.

COTTONOID® Material

COTTONOID compressed rayon provides the natural softness, absorbency and pliability of cotton. A special process interweaves individual rayon fibres to resist shedding, stretching or linting – whether wet or dry.

Soft and Gentle

SUTUREWELD®'s strong bond holds the indicator string to the patty so the surface which lies next to the tissue is completely smooth and free of stitches.

Pliable and Porous



COTTONOID material offers the advantage of being porous so CODMAN Patties and Strips easily conform to the intricate contours of delicate tissue. Yet, you can safely suction through them without aspirating the tissue or the patty strip.

X-Ray Detectable

Marking is where it should be – on the patty or strip.

Unique Packaging

A transparent, sterile envelope shows the size and shape of patties and strips, and the Count Verification Card instantly shows how many patties you started with and how many were used.



Sponge, BICOL

80-1460 BICOL* Collagen Sponge, sterile.
The BICOL Sponges are 3" (76mm) x 4" (102mm) and packaged individually sterile, sold ten packages per box.



Strips, Neurosurgical

CODMAN Surgical Strips, sterile, X-ray detectable, 10 strips per package, 20 packages per dispenser box.

80-1449	1/8" x 6"	(3 x 152mm)
80-1450	1/4" x 6"	(6 x 152mm)
80-1451	1/2" x 6"	(13 x 152mm)
80-1452	3/4" x 6"	(19 x 152mm)
80-1453	1" x 6"	(25 x 152mm)
80-1454	1 1/2" x 6"	(38 x 152mm)
80-1455	2" x 6"	(51 x 152mm)
80-1456	3" x 6"	(76 x 152mm)
80-1457	3 1/2" x 6"	(89 x 152mm)

