

PDS[®] II (Polydioxanone) Monofilament

SUTURE

Description: **Monofilament** For smooth passage through tissue.
Synthetic For minimal tissue reaction.
Absorbable By hydrolysis. Essentially complete in 182-238 days.

Tensile strength: Post implantation Approximate % original strength remaining

14 days	80% (3-0 and larger) 60% (4-0 and smaller)
28 days	70% (3-0 and larger) 40% (4-0 and smaller)
42 days	60% (3-0 and larger) 35% (4-0 and smaller)

Color: Violet or undyed.


















Range: 7-0 to 2 (USP). Supplied as needled sutures and ligatures.

Indications: PDS[®] II sutures are intended for use in general soft tissue approximation, including use in pediatric cardiovascular tissue, in microsurgery and in ophthalmic surgery. PDS[®] II sutures are particularly useful where the combination of an absorbable suture and extended wound support is desirable. Typical areas of use include:











































- Abdominal wall closure
- Intestinal anastomoses
- Ligament and tendon repair



Full details in the Instruction for Use included in every package.












1/2 Circle Taper Point Needles							
Needle	Gauge	Suture Color	Length (cm)	Other information	Number per box	Code	
 TF-6 TF-6	6,5 mm	7-0	●	70		24	Z1370E
 TF-1 TF-1	10 mm	6-0	●	70		24	PUU2971E
 RB-2 RB-2	13 mm	5-0	●	70		36	Z987H
 TF	13 mm	6-0	●	70		36	Z926H
		5-0	●	70		36	Z925H
		4-0	●	70		36	Z924H
 TF TF	13 mm	5-0	●	70		36	Z935H
 RB-1	17 mm	5-0	●	70		36	Z303H
		4-0	●	70		36	Z304H
		3-0	●	70		36	Z305H
 RB-1 RB-1	17 mm	4-0	●	70		36	Z864H
 SH-1	22 mm	4-0	●	70		36	Z310H
		3-0	●	70		36	Z311H
 SH-1 SH-1	22 mm	4-0	●	70		36	Z1310H
		3-0	●	70		36	Z1311H
 JB-1 VISI-BLACK	22 mm	4-0	●	70	VISI-BLACK*	36	Z3100H
		3-0	●	70	VISI-BLACK*	36	Z3110H
 JB-1 JB-1 VISI-BLACK	22 mm	4-0	●	75	VISI-BLACK*	36	Z1312H
 SH	26 mm	4-0	●	70		36	Z315H
		3-0	●	70		36	Z316H
		2-0	●	70		36	Z317H
 SH SH	26 mm	3-0	●	70		36	Z1316H
		2-0	●	90		24	PPF2994E
 JB VISI-BLACK	26 mm	3-0	●	70	VISI-BLACK*	36	Z3160H
		2-0	●	70	VISI-BLACK*	36	Z3170H
 JB JB VISI-BLACK	26 mm	4-0	●	90	VISI-BLACK*	24	X976E
		3-0	●	90	VISI-BLACK*	36	Z5270H
 MH-1	31 mm	2-0	●	150	Loop	36	PV2924SH
 MH	36 mm	3-0	●	70		36	Z322H
		2-0	●	70		36	Z323H

PDS* II
(Polydioxanone) Monofilament


 1/2 Circle Taper Point Needles (Heavier Wire Diameter)						
Needle	Gauge	Suture Color	Length (cm)	Other information	Number per box	Code
 UCL 13 mm	3-0		70		36	Z113H
	2-0		70		36	Z114H
	0		70		36	Z115H
 UCLX 16 mm	3-0		70		36	Z13H
	2-0		70		36	Z14H
	0		70		36	Z15H
 CT-3 22 mm	2-0		70		36	Z328H
	0		70		36	Z329H
 CT-2 26 mm	2-0		70		36	Z333H
	0		70		36	Z334H
	1		70		36	Z335H
 CT-1 36 mm	3-0		70		36	Z338H
	2-0		70		36	Z339H
	2-0		150		12	Z1927G
	0		70		36	Z340H
	1		70		36	Z341H
 MO-2 40 mm	0		150	Loop	24	W9236T
	1		150	Loop	24	W9237T
 CT 40 mm	2-0		90		24	Z357E
	0		90		24	Z358E
	1		90		24	Z359E
 CTX 48 mm	0		90		24	Z370E
	0		150		12	Z1928G
	0		150	Loop	12	PN2925SG
	0		150	Loop	24	W9261T
	1		90		24	Z371E
	1		150	Loop	24	W9262T
	0		150	Loop	12	PS2925SG
 TP-1 65 mm	0		150	Loop	12	PS2925SG
	1		150		12	PS2926G
 XLH 70 mm	0		70		12	Z582G




1/2 Circle ETHIGUARD* Blunt Point Needles							
Needle	Gauge	Suture Color	Length (cm)	Other information	Number per box	Code	
 BT-3	45 mm	1	●	90		24	W9985T
 CTXB	48 mm	0	●	150	Loop	24	W9986T
		1	●	150	Loop	24	W9987T






3/8 Circle Taper Point Needles							
Needle	Gauge	Suture Color	Length (cm)	Other information	Number per box	Code	
 BV175-8 VISI-BLACK	9,3 mm	7-0	●	70	VISI-BLACK*	24	Z161E
 BV-1	9,3 mm	7-0	●	70		24	Z1701E
		6-0	●	70		36	Z1702H
 BV	11 mm	6-0	●	70		36	Z1002H
 C-1	13 mm	6-0	●	70		36	Z1012H
		5-0	●	70		36	Z1013H
 C-1	13 mm	6-0	●	70		36	Z1721H
		5-0	●	70		36	Z1722H
 BB	17 mm	5-0	●	90		24	KCC2992E


5/8 Circle Taper Point Needles							
Needle	Gauge	Suture Color	Length (cm)	Other information	Number per box	Code	
 UR-6	26 mm	2-0	●	70		36	VC2973H
 UR-6	26 mm	3-0	●	70		24	Z805E
 MF-110 HVY	110 mm	0	●	150	Loop	12	W9168









PDS* II
(Polydioxanone) Monofilament

5/8 Circle ETHIGUARD* Blunt Point Needles						
Needle	Gauge	Suture Color	Length (cm)	Other information	Number per box	Code
 <p>UB-110 110 mm</p>	0	●	240	Loop	12	W9968




Straight Taper Point Needles						
Needle	Gauge	Suture Color	Length (cm)	Other information	Number per box	Code
 <p>ST-4 19 mm</p>	5-0	●	35		36	BA2932H
 <p>ST-4 19 mm</p>	4-0	●	35		24	Z620E
 <p>STA 65 mm</p>	0	●	70		24	Z630E

1/2 Circle TAPERCUT* Needles						
Needle	Gauge	Suture Color	Length (cm)	Other information	Number per box	Code
 <p>V-5 17 mm</p>	5-0	●	70		36	Z995H
	4-0	●	70		36	Z994H
	3-0	●	70		36	Z276H
	2-0	●	70		36	Z274H
 <p>V-6 22 mm</p>	4-0	●	70		36	Z275H
	3-0	●	70		36	Z277H
 <p>V-7 26 mm</p>	4-0	●	70		36	Z997H
	3-0	●	70		36	Z998H
	2-0	●	70		36	Z999H
 <p>V-37 40 mm</p>	0	●	150	Loop	24	Z1932
 <p>V-40 48 mm</p>	0	●	150	Loop	24	Z1933E



 **3/8 Circle CC-Needles**

Needle	Gauge	Suture Color	Length (cm)	Other information	Number per box	Code
 CC 9,3 mm	7-0		70		24	Z1751E
 CC-CC 9,3 mm	7-0		70		24	Z1711E
	6-0		70		36	Z1712H
 CC-1-CC-1 13 mm	6-0		70		36	Z1032H
	5-0		70		24	Z1033E



 **3/8 Circle TAPERCUT* Needles**

Needle	Gauge	Suture Color	Length (cm)	Other information	Number per box	Code
 V-18 13 mm	6-0		70		36	Z991H
	5-0		70		36	Z990H







 **Fish hook TAPERCUT* Needles**





Needle	Gauge	Suture Color	Length (cm)	Other information	Number per box	Code
 V-38 31 mm	1		70		24	W9334T




 **Straight TAPERCUT* Needles**




Needle	Gauge	Suture Color	Length (cm)	Other information	Number per box	Code
 V-47 38 mm	4-0		35		24	Z425E


PDS* II
(Polydioxanone) Monofilament


▼ 1/2 Circle Reverse Cutting Needles (Heavier Wire Diameter)							
Needle	Gauge	Suture Color	Length (cm)	Other information	Number per box	Code	
 OS-4	22 mm	1	●	70	Heavy Needle	24	W9286T
 X-1	22 mm	3-0	○	70		24	Z458E
 CP-2	26 mm	2-0	●	70		36	Z969H
		0	●	70		36	Z983H
		1	●	70		36	Z984H
 CP-1	36 mm	2-0	●	70		36	Z466H
		0	●	70		36	Z467H
		1	●	70		36	Z468H
 CP	40 mm	0	●	90		24	Z485E
		1	●	90		24	W9324T
		2	●	70		24	MH2976E
 CPX	48 mm	0	●	90		24	Z1057E
		1	●	90		24	W9299T

▼ 3/8 Circle Reverse Cutting PRIME* Needles							
Needle	Gauge	Suture Color	Length (cm)	Other information	Number per box	Code	
 P-1	11 mm	6-0	●	45		36	Z487H
		6-0	○	45		24	Z489E
		6-0	○	45	MULTIPASS*	36	MPZ489H
		5-0	○	45		24	Z490E
		5-0	○	45	MULTIPASS*	36	MPZ490H
 P-3	13 mm	5-0	●	45		36	Z463H
		5-0	○	45		24	Z493E
		5-0	○	45	MULTIPASS*	36	MPZ493H
		4-0	●	45		36	Z464H
		4-0	○	45		24	Z494E
		4-0	○	45	MULTIPASS*	36	MPZ494H
 PS-3	16 mm	5-0	○	45	MULTIPASS*	36	MPZ500H
		5-0	○	45		24	W9872T
		4-0	○	45	MULTIPASS*	36	MPZ501H
		4-0	○	45		24	W9873T
 PS-2	19 mm	4-0	○	45		24	W9874T
		4-0	○	45	MULTIPASS*	36	MPZ496H
		3-0	○	45	MULTIPASS*	36	MPZ497H
		3-0	○	45		24	Z497E



▼ 3/8 Circle Reverse Cutting PRIME* Needles							
Needle	Gauge	Suture Color	Length (cm)	Other information	Number per box	Code	
 PS-2 S	19 mm	4-0	○	45	Slimmed	24	Z498E
		4-0	●	45	Slimmed	36	Z513H
 PS-1	24 mm	4-0	○	45	MULTIPASS*	36	MPZ682H
		3-0	○	45	MULTIPASS*	36	MPZ683H
		3-0	○	45		24	Z683E
		2-0	○	45		24	Z684E
 PS	26 mm	2-0	○	45	MULTIPASS*	36	MPZ9626H
		2-0	○	70		24	W9626T

▼ 3/8 Circle Reverse Cutting Needles							
Needle	Gauge	Suture Color	Length (cm)	Other information	Number per box	Code	
 FS-2	19 mm	5-0	○	70		24	Z421E
		4-0	●	30	Loop	36	Z282H
		3-0	○	45		24	Z293E
		3-0	○	70		24	Z423E
 FS-2S	19 mm	4-0	○	45	Slimmed	24	Z292E
		4-0	○	70	Slimmed	24	Z422E
 FS-1	24 mm	4-0	○	70		24	Z441E
		3-0	○	70		24	Z442E
		2-0	○	70		24	Z443E

▲ 3/8 Circle Conventional Cutting PRIME* Needles							
Needle	Gauge	Suture Color	Length (cm)	Other information	Number per box	Code	
 PC-1	13 mm	6-0	○	45		24	W9861T

▼ 1/2 Circle Reverse Cutting PRIME* Needles							
Needle	Gauge	Suture Color	Length (cm)	Other information	Number per box	Code	
 P-2	8 mm	5-0	●	45		36	JG2942H

PDS* II
(Polydioxanone) Monofilament

▼ Straight Cutting Needles							
Needle	Gauge	Suture Color	Length (cm)	Other information	Number per box	Code	
 SC-1 13 mm	4-0	●	45		36	Z998H	
 KS 60 mm	4-0	○	70		36	Z862H	
	3-0	○	70		36	Z863H	
	2-0	○	70		36	Z864H	

Ligatures						
Needle	Gauge	Suture Color	Length (cm)	Other information	Number per box	Code
SUTUPAK*	2-0	●	2x70		36	Z625H
	0	●	2x70		24	Z626E
	1	●	2x70		24	Z627E