

# ETHICON

a Johnson & Johnson company



**7-0**  
(0.5 Ph. Eur.)

**BV175-8**  
9.3 mm 3/8c  
ROUND BODIED



ETHALLOY®  
VISI-BLACK®

**75cm**

**PROLENE\***

Polypropylene  
blu monofilamento sutura non assorbibile  
azul monofilamento sutura no absorbible  
blauw monofil chirurg. hechtmateriaal  
niet resorbeerbaar

36 pezzi / unidades / stuks  
STERILE EO 2

LOT APE048  
2013-07

ATRALOC®

**ETHICON\***

**7-0**  
(0.5 Ph. Eur.)  
BV175-8  
9.3 mm 3/8c  
ROUND BODIED



**ETHICON\***  
**PROLENE\***  
Polypropylene

ETHICON, INC.  
RABBIT CREEK, OHIO, U.S.A.

**PROLENE\***  
Monofilament Polypropylene

# PROLENE\*

Monofilament Polypropylene

## SUTURE

**Description:** **Monofilament** For smooth passage through tissue.  
**Synthetic** For minimal tissue reaction.  
**Non-absorbable** Provides permanent tensile strength retention in tissue.

**Color:** Blue or clear.

**Range:** 10-0 to 1 (USP). Supplied as needed sutures, armed with, for example, VISI-BLACK\* needles, MULTIPASS\* needles or CC-needles. Also available with pledgets.

PROLENE\* sutures are also available as HEMO-SEAL\*, with suture diameter to needle diameter ratio close to 1. HEMO-SEAL\* is particularly useful in reducing needle hole leakage (for example in synthetic vascular grafts).

**Indications:** PROLENE\* sutures are intended for use in general soft tissue approximation and/or ligation, including use in cardiovascular, ophthalmic and neurosurgical procedures. Typical areas of use include:

- Cardiovascular anastomoses
- Subcuticular sutures

Full details in the Instruction for Use included in every package.



# PROLENE\* Monofilament Polypropylene

## 3/8 Circle Taper Point Needles

Needle		Gauge	Suture Color	Length (cm)	Other information	Number per box	Code
BV130-5 	6,5 mm	8-0	BLUE	7	VISI-BLACK*	24	EH7456E
		8-0	BLUE	45	EVERPOINT*	36	EP8730H
		8-0	BLUE	45	VISI-BLACK*	12	X916G
		8-0	BLUE	45	MULTIPASS*	24	EH7470E
		8-0	BLUE	45	MULTIPASS*	12	W2777
		8-0	BLUE	60	VISI-BLACK*	12	F1894
		8-0	BLUE	60	EVERPOINT*	36	EP8732H
		8-0	BLUE	60	MULTIPASS*	36	KHH5662H
		8-0	BLUE	75	VISI-BLACK*	24	EH7459E
		8-0	BLUE	75	MULTIPASS*	24	EH7469E
		8-0	BLUE	75	Straight pack, MULTIPASS*	12	EH7469LG
		7-0	BLUE	60	VISI-BLACK*	36	8207H
		7-0	BLUE	2x75	Straight pack, MULTIPASS*	12	KHL5673JG
BV130-5 	6,5 mm	8-0	BLUE	13		12	W2775
BV175-6 	8 mm	8-0	BLUE	45	MULTIPASS*	36	8740H
		8-0	BLUE	60	VISI-BLACK*	36	8841H
		8-0	BLUE	60	EVERPOINT*	36	EP8741H
		8-0	BLUE	60	MULTIPASS*	36	8741H
		7-0	BLUE	60	VISI-BLACK*	36	8335H
		7-0	BLUE	60	MULTIPASS*	36	8735H
		7-0	BLUE	60	EVERPOINT*	36	EP8735H
		7-0	BLUE	2X60	EVERPOINT*, Straight pack	12	EPM8737
		7-0	BLUE	2x60	Straight pack, VISI-BLACK*	12	M8935
		7-0	BLUE	2x60	Straight pack, MULTIPASS*	12	M8737G
BV175-7 	8 mm	8-0	BLUE	2x60	Straight pack, VISI-BLACK*	12	M8967
		7-0	BLUE	60		36	8766H
		7-0	BLUE	2x60	Straight pack	12	M8766
BV175-8 	9,3 mm	7-0	BLUE	45	VISI-BLACK*	36	EH7225H
		7-0	BLUE	45	EVERPOINT*	36	EP8746H
		7-0	BLUE	45	MULTIPASS*	36	8701H
		7-0	BLUE	60	EVERPOINT*	36	EP8747H
		7-0	BLUE	60	MULTIPASS*	12	W8702
		7-0	BLUE	4x60	VISI-BLACK*	12	X967G
		7-0	BLUE	4X60	EVERPOINT*, Straight pack	12	EPM8745
		7-0	BLUE	4x60	Straight pack, MULTIPASS*	12	W8702S
		7-0	BLUE	75	VISI-BLACK*	36	EH7222H
		7-0	BLUE	75	EVERPOINT*	36	EP8755H
		7-0	BLUE	75	MULTIPASS*	36	EH7405H
		7-0	BLUE	2x75	Straight pack, VISI-BLACK*	12	MEH7224JG
		7-0	BLUE	4x75	Straight pack, MULTIPASS*	12	M8702G






# PROLENE\* Monofilament Polypropylene

## 3/8 Circle Taper Point Needles






Needle		Gauge	Suture Color	Length (cm)	Other information	Number per box	Code
	9,3 mm	7-0	BLUE	75	MULTIPASS*	36	EH7402H
	9,3 mm	7-0	BLUE	60	VISI-BLACK*	12	W8304
		7-0	BLUE	60	EVERPOINT*	36	EP8702H
		7-0	BLUE	2x60	Straight pack, VISI-BLACK*	12	KAA5663LG
		6-0	BLUE	60	VISI-BLACK*	12	W8305
		6-0	BLUE	60	EVERPOINT*	36	EP8805H
		6-0	BLUE	75	VISI-BLACK*	36	EH7223H
		6-0	BLUE	75	EVERPOINT*	36	EP8709H
		6-0	BLUE	75		36	EH7406H
5-0	BLUE	75	EVERPOINT*, HEMO-SEAL*	36	EPH9702H		
	9,3 mm	6-0	BLUE	75		36	EH7403H
	9,3 mm	7-0	BLUE	75	VISI-BLACK*, Heavier Wire Diameter	36	EH7238H
		6-0	BLUE	75	VISI-BLACK*, Heavier Wire Diameter	36	EH7240H
		6-0	BLUE	75	Heavier Wire Diameter	36	EH7245H
	11 mm	6-0	BLUE	75		36	8776H
	13 mm	7-0	BLUE	60		12	W8725
		7-0	BLUE	75		36	8708H
		7-0	BLUE	75	EVERPOINT*	36	EP8708H
		6-0	BLUE	45	VISI-BLACK*	36	X577H
		6-0	BLUE	60	VISI-BLACK*	12	W8307
		6-0	BLUE	60		12	W8718
		6-0	BLUE	75	VISI-BLACK*	36	EH7228H
		6-0	BLUE	75		36	8706H
		6-0	BLUE	75	EVERPOINT*	36	EP8706H
		6-0	BLUE	4x75	EVERPOINT*, Straight pack	12	EPM8706
		6-0	BLUE	4x75	Straight pack, VISI-BLACK*	12	M807C
		6-0	BLUE	90	VISI-BLACK*	36	EH7230H
		5-0	BLUE	45		36	8717H
		5-0	BLUE	60		36	KBB5661H
		5-0	BLUE	60	EVERPOINT*, HEMO-SEAL*	36	EPH8725H
		5-0	BLUE	70	With pledget 4x3mm	12	KBB5671VG
		5-0	BLUE	75		36	EH7477H
		5-0	BLUE	75	EVERPOINT*, HEMO-SEAL*	36	EPH7477H
		5-0	BLUE	75	VISI-BLACK*	36	EH7229H
		5-0	BLUE	2x75	Straight pack, VISI-BLACK*	12	MEH7233JG
		5-0	BLUE	90	VISI-BLACK*	36	EH7231H
		5-0	BLUE	90		36	8720H
		5-0	BLUE	2x90	Straight pack	12	KBL5691JG
		5-0	BLUE	90	EVERPOINT*, HEMO-SEAL*	36	EPH8720H
4-0	BLUE	75		36	EH7478H		

# PROLENE\* Monofilament Polypropylene

## 3/8 Circle Taper Point Needles










Needle		Gauge	Suture Color	Length (cm)	Other information	Number per box	Code
	13 mm	6-0	BLUE	75	VISI-BLACK*	36	EH7226H
		5-0	BLUE	75	VISI-BLACK*	36	EH7227H
		6-0	BLUE	75	EVERPOINT*	36	EP8889H
		6-0	BLUE	75		36	8889H
		5-0	BLUE	75		36	8890H
		5-0	BLUE	75	EVERPOINT*, HEMO-SEAL*	36	EPH8890H
		4-0	BLUE	75		36	EH7476H
	13 mm	5-0	BLUE	75	VISI-BLACK*, Heavier Wire Diameter	36	EH7257H
		5-0	BLUE	75	Heavier Wire Diameter	36	EH7255H
	13 mm	5-0	BLUE	75	Heavier Wire Diameter	36	EH7254H
	17 mm	5-0	BLUE	90		36	8580H
		4-0	BLUE	90		36	8581H
	17 mm	5-0	BLUE	75		36	8880H
		4-0	BLUE	75		36	8881H

## 1/2 Circle Taper Point Needles

Needle		Gauge	Suture Color	Length (cm)	Other information	Number per box	Code
	6,5 mm	7-0	BLUE	45		24	EH7813E
		7-0	BLUE	75		36	PQQ5671H
		6-0	BLUE	45		24	EH7814E
		6-0	BLUE	60		24	PQQ5660E
	10 mm	6-0	BLUE	75		36	EH7835H
		5-0	BLUE	60	HEMO-SEAL*	36	PUU5662SH
	11 mm	7-0	BLUE	60		36	8206H
	13 mm	6-0	BLUE	45		36	8714H
		6-0	BLUE	75	EVERPOINT*	36	EP8711H
		6-0	BLUE	75		36	8711H
		5-0	BLUE	75		36	8710H
		5-0	BLUE	4x75	Straight pack	12	M8710G
		5-0	BLUE	75	EVERPOINT*, HEMO-SEAL*	36	EPH8710H
		5-0	BLUE	90	EVERPOINT*, HEMO-SEAL*	36	EPH8716H
	13 mm	6-0	BLUE	75		36	EH7400H
		5-0	BLUE	60		36	8205H
		5-0	BLUE	75		36	EH7401H
		5-0	BLUE	75	With 2x2 mm firm pledget	36	EH7412H
		4-0	BLUE	90		36	EH7411H

# PROLENE\* Monofilament Polypropylene

## 1/2 Circle Taper Point Needles

Needle		Gauge	Suture Color	Length (cm)	Other information	Number per box	Code
	17 mm	5-0	BLUE	75	HEMO-SEAL*	36	HS6856H
		5-0	BLUE	90	VISI-BLACK*	36	VB50H
		5-0	BLUE	90		36	8556H
		5-0	BLUE	90	HEMO-SEAL*	36	HS8556H
		4-0	BLUE	90	VISI-BLACK*	36	8357H
		4-0	BLUE	90		36	8557H
		4-0	BLUE	90	HEMO-SEAL*	36	HS6857H
		3-0	BLUE	90		36	8558H
	17 mm	5-0	BLUE	75		36	8870H
		4-0	BLUE	75		36	8871H
		3-0	BLUE	75		36	8872H
	20 mm	4-0	BLUE	90	HEMO-SEAL*	36	HS6861H
	22 mm	4-0	BLUE	90		36	EH7585H
		4-0	BLUE	90	With pledget 5x8mm soft	36	PEE5693SH
		3-0	BLUE	90		36	EH7586H
	22 mm	4-0	BLUE	75		36	EH7581H
		3-0	BLUE	75		36	EH7582H
	26 mm	4-0	BLUE	90	VISI-BLACK*	36	EH8087H
		4-0	BLUE	90	HEMO-SEAL*	36	HS6855H
		4-0	BLUE	90	Pledget 6x3mm firm, Straight pack	12	EH7190LG
		4-0	BLUE	90		36	8521H
		4-0	BLUE	90	Pledget 5x8mm soft, VISI-BLACK*	36	EH8062H
		4-0	BLUE	120		36	EH7583H
		3-0	BLUE	90	VISI-BLACK*	36	EH8086H
		3-0	BLUE	90		36	8522H
		3-0	BLUE	90	HEMO-SEAL*	36	HS6822H
		3-0	BLUE	120		36	EH7584H
		2-0	BLUE	90		36	8523H
	26 mm	4-0	BLUE	75		36	8831H
		3-0	BLUE	75		36	8832H
		2-0	BLUE	75		36	8833H
		0	BLUE	75		36	8834H
	31 mm	3-0	BLUE	120	VISI-BLACK*	36	X1295H
		3-0	BLUE	90	VISI-BLACK*	36	EH8091H
	36 mm	3-0	BLUE	90		36	EH7499H
		2-0	BLUE	90		36	PGG5694H
		2-0	BLUE	122		36	8853H

# PROLENE\* Monofilament Polypropylene


## 1/2 Circle Taper Point Needles (Heavier Wire Diameter)

Needle		Gauge	Suture Color	Length (cm)	Other information	Number per box	Code
 CT-2	26 mm	2-0	BLUE	75		36	8411H
		0	BLUE	75		36	8412H
		1	BLUE	90		36	PK5696H
 CT-1	36 mm	3-0	BLUE	75		36	8422H
		2-0	BLUE	75		36	8423H
		0	BLUE	75		36	8424H
		1	BLUE	75		36	8425H
 CT	40 mm	0	BLUE	75		36	8434H
		1	BLUE	100		12	W742
 CTX	48 mm	0	BLUE	90		36	PN5695H
		1	BLUE	150	Loop	12	8527G

## Compound Curve Taper Point Needles

Needle		Gauge	Suture Color	Length (cm)	Other information	Number per box	Code
 ABV-1	 ABV-1	9,3 mm	7-0	BLUE	75	36	EH8000H
 AC-1	 AC-1	13 mm	6-0	BLUE	75	36	EH8010H
			5-0	BLUE	75	36	EH8014H

## 5/8 Circle Taper Point Needles








Needle		Gauge	Suture Color	Length (cm)	Other information	Number per box	Code
 UR-6	26 mm	2-0	BLUE	75		36	VC5673H

## Straight Taper Point Needles

Needle		Gauge	Suture Color	Length (cm)	Other information	Number per box	Code
 ST-4	19 mm	4-0	BLUE	35		36	EH7508H
 ST-1	45 mm	3-0	BLUE	35		36	EH7921H
 ST-70	70 mm	2-0	BLUE	75		12	W8400

# PROLENE\* Monofilament Polypropylene








## 3/8 Circle CC-Needles

Needle		Gauge	Suture Color	Length (cm)	Other information	Number per box	Code
CC175-6 	8 mm	8-0	BLUE	60	VISI-BLACK*	12	W8101
CC175-8 	9,3 mm	8-0	BLUE	60		24	EH7972E
		8-0	BLUE	60	MULTIPASS*	12	W8703
		8-0	BLUE	60	EVERPOINT*	36	EP7972H
		7-0	BLUE	45		36	8705H
		7-0	BLUE	75	EVERPOINT*	36	EP8812SLH
		7-0	BLUE	75	MULTIPASS*	36	8704H
6-0	BLUE	60	MULTIPASS*	12	W8807		
CC175-8 	9,3 mm	7-0	BLUE	75		36	EH7970H
CC 	9,3 mm	7-0	BLUE	75	VISI-BLACK*	36	EH7237H
		7-0	BLUE	75		36	EH7236H
		7-0	BLUE	75	EVERPOINT*	36	EP8812H
		6-0	BLUE	75	VISI-BLACK*	36	EH7234H
		6-0	BLUE	75	EVERPOINT*	36	EP8811H
		6-0	BLUE	75		36	EH7235H
		6-0	BLUE	60	EVERPOINT*	36	EP8807H
		5-0	BLUE	60	HEMOSEAL*	12	W8662
CC-11 	11 mm	6-0	BLUE	60		12	W8802
CC-1 	13 mm	6-0	BLUE	45		36	8719H
		6-0	BLUE	75	MULTIPASS*	36	8707H
		6-0	BLUE	2x75	MULTIPASS*, Straight pack	12	X535G
		5-0	BLUE	75		36	EH7977H
		5-0	BLUE	75		36	EH7251H
		5-0	BLUE	75	MULTIPASS	36	7475H
		5-0	BLUE	90		36	8721H
CC-1 	13 mm	6-0	BLUE	75		36	EH7975H
		5-0	BLUE	75		36	EH7976H








# PROLENE\* Monofilament Polypropylene

## 3/8 Circle TAPERCUT\* Needle

Needle	Gauge	Suture Color	Length (cm)	Other information	Number per box	Code	
 MV-14	5,5 mm	8-0	BLUE	60		24	EH7828E
 V-18	13 mm	6-0	BLUE	75		36	EH7474H
		5-0	BLUE	75		36	EH7475H
 V-18	13 mm	6-0	BLUE	75		36	EH7471H
		5-0	BLUE	75		36	EH7472H
 V-18 X	13 mm	5-0	BLUE	75	Heavier Wire Diameter	36	EH7253H
 V-4	17 mm	5-0	BLUE	75		36	F2829
		4-0	BLUE	75		36	F2828
		4-0	BLUE	90		36	8925H
 V-4	17 mm	5-0	BLUE	75		36	8941H
		4-0	BLUE	75		36	8942H
		3-0	BLUE	75		36	8943H
 V-26	26 mm	3-0	BLUE	90		36	EH7788H

## 1/2 Circle CC-Needles

Needle	Gauge	Suture Color	Length (cm)	Other information	Number per box	Code	
 CC-6	11 mm	7-0	BLUE	60	VISI-BLACK*	12	X1004G
 CC-4	13 mm	6-0	BLUE	60	MULTIPASS*	12	W8815
		6-0	BLUE	75	VISI-BLACK*	36	EH7242H
		5-0	BLUE	60	HEMO-SEAL*	12	W8664
		5-0	BLUE	75	VISI-BLACK*	36	EH7241H
		5-0	BLUE	2x75	VISI-BLACK*, Straight pack	12	MEH7241JG
		5-0	BLUE	75		12	W8816
		4-0	BLUE	75	HEMO-SEAL*	12	W8665
 CC-16	16 mm	3-0	BLUE	90	HEMO-SEAL*	12	W8667
 CC-25	26 mm	4-0	BLUE	90		12	W8845
 CC-30	31 mm	3-0	BLUE	90		12	W8849

# PROLENE\* Monofilament Polypropylene

## Compound Curve CC-Needles


Needle		Gauge	Suture Color	Length (cm)	Other information	Number per box	Code
	9,3 mm	7-0	BLUE	60	VISI-BLACK*	36	EH8020H
		7-0	BLUE	75	VISI-BLACK*	24	EH8021E
		6-0	BLUE	75	DOUBLE ARMED ACC + ACC-1	36	EH8035H
		6-0	BLUE	75	VISI-BLACK*	36	EH8025H
	13 mm	6-0	BLUE	75	VISI-BLACK*	36	EH8030H
		6-0	BLUE	90	VISI-BLACK*	36	EH8031H
		6-0	BLUE	2x75	VISI-BLACK*, Straight pack	12	MEH8034JG
		5-0	BLUE	75	VISI-BLACK*	36	EH8036H
		5-0	BLUE	90	VISI-BLACK*	36	EH8037H
	13 mm	6-0	BLUE	75	VISI-BLACK*	36	EH8032H

## 1/2 Circle TAPERCUT\* Needles










Needle		Gauge	Suture Color	Length (cm)	Other information	Number per box	Code
	17 mm	5-0	BLUE	90		36	8934H
		4-0	BLUE	90		36	8935H
		4-0	BLUE	90	Straight pack, Pledgets	12	8954G
		3-0	BLUE	90		36	8936H
		2-0	BLUE	90		36	8937H
	17 mm	5-0	BLUE	75		36	8951H
		4-0	BLUE	75		36	8952H
		3-0	BLUE	75		36	8953H
	22 mm	4-0	BLUE	90		36	EH7774H
		3-0	BLUE	90		36	EH7775H
	22 mm	2-0	BLUE	75		36	EH7772H
	26 mm	4-0	BLUE	90		36	8975H
		3-0	BLUE	90		36	8976H
		3-0	BLUE	90	Pledgets	36	OBB5694SH
		3-0	BLUE	120		36	EH7778H
		2-0	BLUE	90		36	8977H
		2-0	BLUE	4x90	Straight pack, Pledgets	12	OBL5595LG
		2-0	BLUE	4x90	Straight pack	12	OBL5695LG
		2-0	BLUE	120		36	EH7779H
	26 mm	4-0	BLUE	75		36	8962H
		3-0	BLUE	75		36	8963H
		2-0	BLUE	75		36	8964H
	36 mm	0	BLUE	75		36	8444H

# PROLENE\* Monofilament Polypropylene

## Straight TAPERCUT\* Needles








Needle	Gauge	Suture Color	Length (cm)	Other information	Number per box	Code
 V-47	38 mm	4-0	BLUE	35		36 8592H

## 3/8 Circle Reverse Cutting PRIME\* Needles




Needle	Gauge	Suture Color	Length (cm)	Other information	Number per box	Code
 P-6	8 mm	7-0	BLUE	45		36 8648H
		6-0	BLUE	45		24 W8870T
 P-1	11 mm	7-0	BLUE	75	MULTIPASS*	36 MPP2854H
		6-0	BLUE	45		36 8697H
		6-0	BLUE	45	MULTIPASS*	36 MPP8697H
 P-3	13 mm	6-0	BLUE	45		36 8695H
		6-0	BLUE	45	MULTIPASS*	36 MPP8695H
		5-0	BLUE	45	MULTIPASS*	36 8698H
		5-0	UNDYED	75	MULTIPASS*	36 MPP2858H
		4-0	BLUE	45	MULTIPASS*	36 8699H
		4-0	UNDYED	45		36 JA5742H
		4-0	UNDYED	45	MULTIPASS*	36 MPP8604H
 PS-3	16 mm	6-0	BLUE	75	MULTIPASS*	36 MPP2832H
		5-0	BLUE	45		36 8681H
		5-0	BLUE	45	MULTIPASS*	36 MPP8681H
		4-0	BLUE	45		36 EH7764H
		4-0	BLUE	45	MULTIPASS*	36 MPP7764H
 PS-2		5-0	BLUE	45	MULTIPASS*	36 8686H
		4-0	BLUE	45		24 W8015T
		4-0	BLUE	45	MULTIPASS*	36 MPP8682H
		4-0	BLUE	75	MULTIPASS*	36 EH7167H
		4-0	UNDYED	75	MULTIPASS*	36 MPP2861H
		3-0	BLUE	60	MULTIPASS*	36 EH7691H
		3-0	UNDYED	75	MULTIPASS*	36 MPP2862H
 PS-2 S	19 mm	4-0	BLUE	45	Slimmed	36 8682H
 PS-1	24 mm	4-0	BLUE	45		36 EH7758H
		4-0	BLUE	75	MULTIPASS*	36 MPP7695H
		3-0	BLUE	45	MULTIPASS*	12 8663G
		3-0	BLUE	60		12 EH7666G
		3-0	BLUE	75		36 EH7568H
 PS	26 mm	2-0	BLUE	45	MULTIPASS*	36 MPP8026H
 PSL	30 mm	2-0	BLUE	75		36 EH7121H

# PROLENE\* Monofilament Polypropylene


## 3/8 Circle Reverse Cutting Needles

Needle	Gauge	Suture Color	Length (cm)	Other information	Number per box	Code	
 C-2	13 mm	6-0	BLUE	45	250 microns	36	FD5640H
 FS-3	16 mm	6-0	BLUE	45		36	8660H
		5-0	BLUE	45		36	8614H
 FS-2	19 mm	5-0	BLUE	45		36	8661H
		4-0	BLUE	45		36	EH7151H
		4-0	BLUE	75		36	EH7150H
		3-0	BLUE	45		36	8665H
		3-0	BLUE	75		36	EH7693H
 FS-1	24 mm	4-0	BLUE	45		36	8629H
		3-0	BLUE	45		36	8684H
		3-0	BLUE	75		36	EH7694H
 FS	26 mm	2-0	BLUE	45		36	8685H
 FSL	30 mm	3-0	BLUE	75		36	8675H
		2-0	BLUE	75		36	EH7697H
		0	BLUE	45		36	EH7920H
 FSLX	36 mm	3-0	BLUE	75		36	8649H
		2-0	BLUE	75		36	8689H
		0	BLUE	75		36	EH7561H

## 3/8 Circle Conventional Cutting PRIME\* Needles

Needle	Gauge	Suture Color	Length (cm)	Other information	Number per box	Code	
 PC-3	16 mm	6-0	BLUE	45		24	W8873T
		6-0	BLUE	45	MULTIPASS*	36	MPP8873H
		5-0	BLUE	45	MULTIPASS*	36	8635H
		4-0	BLUE	45		12	W525
		4-0	BLUE	45		24	W8007T
		4-0	BLUE	45	MULTIPASS*	36	MPP8875H
 PC-5	19 mm	5-0	BLUE	45	MULTIPASS*	36	MPP8876H
		4-0	BLUE	45		24	W8877T
		4-0	BLUE	45	MULTIPASS*	36	MPP8877H
 PC-25	26 mm	5-0	BLUE	45		24	W8881T
		4-0	BLUE	45		24	W8879T
		3-0	BLUE	45		12	W538

## 1/2 Circle Reverse Cutting PRIME\* Needles

Needle	Gauge	Suture Color	Length (cm)	Other information	Number per box	Code	
 PS-4	16 mm	4-0	BLUE	45		36	EH7759H

# PROLENE\* Monofilament Polypropylene

## 1/2 Circle Reverse Cutting Needles

Needle	Gauge	Suture Color	Length (cm)	Other information	Number per box	Code
 J-1	17 mm	3-0	BLUE	45		36 EH7545H
 CP-2	26 mm	2-0	BLUE	75		36 F2819



## Straight Cutting Needles

Needle	Gauge	Suture Color	Length (cm)	Other information	Number per box	Code
 SC-1	13 mm	4-0	BLUE	45		36 EH7730H
 KS	60 mm	3-0	BLUE	75		36 8622H
 KS		2-0	BLUE	75		36 8623H

## 3/8 Circle Taper Point Needles Special





Needle	Gauge	Suture Color	Length (cm)	Other information	Number per box	Code
 BV75-3	3,8 mm	10-0	BLUE	13		12 W2790
 BV130-4	4,7 mm	9-0	BLUE	13		12 W2780
		8-0	BLUE	13		12 W2779
 BV-6	4,8 mm	10-0	BLUE	10		12 EH7595G
 BV-4	4,9 mm	10-0	BLUE	15		12 EH7596G
 BV100-4	5,1 mm	10-0	BLUE	13		12 W2794
		9-0	BLUE	13		12 W2783
 BV100-4  BV100-4	5,1 mm	9-0	BLUE	45		12 F1891
 BV-2	6,5 mm	9-0	BLUE	15		12 EH7597G

## 1/4 Circle TAPERCUT\* Needles


Needle	Gauge	Suture Color	Length (cm)	Other information	Number per box	Code
 CIF-4  CIF-4	13 mm	10-0	BLUE	20		12 788G

# PROLENE\* Monofilament Polypropylene


## 3/8 Circle MICROPOINT\* Spatula

Needle	Gauge	Suture Color	Length (cm)	Other information	Number per box	Code	
GS-12 	5,5 mm	10-0	BLUE	30	12	EH7473G	
GS-12 GS-12 	5,5 mm	10-0	BLUE	30	12	EH7598G	
CS-35-C CS-35-C 	6 mm	10-0	BLUE	30	12	EH7984G	
TG140-6 TG140-6 	6,5 mm	10-0	BLUE	30	150 MICRONS	12	W1777
GS-9 GS-9 	6,5 mm	9-0	BLUE	15	12	1754G	



## 1/4 Circle CS ULTIMA\* Spatula

Needle	Gauge	Suture Color	Length (cm)	Other information	Number per box	Code
CTC-6 L 	13 mm	10-0	BLUE	10	12	EH7897G



## Ski MICROPOINT\* Spatula

Needle	Gauge	Suture Color	Length (cm)	Other information	Number per box	Code
CSS-160 CSS-160 	16 mm	10-0	BLUE	30	12	EH7997G

## Straight MICROPOINT\* Spatula

Needle	Gauge	Suture Color	Length (cm)	Other information	Number per box	Code	
STC-6 STC-6 	16 mm	10-0	BLUE	20	12	1713G	
STC-6 	16 mm	10-0	BLUE	15	Loop	12	EH7900G

## Special Combinations

Needle	Gauge	Suture Color	Length (cm)	Other information	Number per box	Code
STC-6 CS-35-C 	16 mm + 6 mm	10-0	BLUE	20	12	ET5631SG
STC-6 GS-9 	6,5 mm + 16 mm	9-0	BLUE	20	12	1753G