

# SILK

Braided Silk / Virgin silk

## SUTURE

### Description:

<b>Braided</b>	For easy handling and secure knot tying.
<b>Coated</b>	With wax for easy handling.
<b>Natural</b>	Produced from cocoons of the silk worm <i>Bombyx Mori</i> .
<b>Non-absorbable</b>	Provides prolonged tensile strength retention in tissue for up to 3 months.

### Color:

Black or white.

### Range:

10-0 to 4 (USP). Supplied as needled sutures or ligatures. Sizes 10-0 to 8-0 are manufactured of virgin silk and are twisted in construction. Virgin silk fibres are held together by natural gum secreted by the silk worm.

### Indications:






SILK sutures are intended for use in general soft tissue approximation and/or ligation, including use in cardiovascular, ophthalmic and neurosurgical procedures. Typical areas of use include:

- Ophthalmic surgery
- Dental surgery
- Ligation





Full details in the Instruction for Use included in every package.




**1/2 Circle Taper Point Needles**

Needle	Gauge	Suture Color	Length (cm)	Other information	Number per box	Code
 TF	13 mm	5-0	●	45	36	N266H
		4-0	●	45	36	N267H
 RB-1	17 mm	5-0	●	75	36	K870H
		4-0	●	75	36	K871H
		3-0	●	75	36	K872H
		2-0	●	75	36	K873H
 SH-1	22 mm	3-0	●	75	36	EH7343H
		2-0	●	75	36	EH7344H
 SH	26 mm	4-0	●	75	36	K831H
		3-0	●	75	36	K832H
		2-0	●	75	36	K833H
		0	●	75	36	K834H
		1	●	75	36	K835H
 MH	36 mm	2-0	●	75	36	K843H
		0	●	75	36	K844H
		1	●	75	36	K845H
		2	●	75	36	K846H

**3/8 Circle Taper Point Needles**



Needle	Gauge	Suture Color	Length (cm)	Other information	Number per box	Code
 BV-1	9,3 mm	7-0	●	30	36	K803H
 C-1	13 mm	6-0	●	75	36	K889H
		5-0	●	75	36	K890H
		4-0	●	75	36	K891H
 C-1 C-1	13 mm	6-0	●	45	36	706H
 BB	17 mm	4-0	●	75	36	K881H


**Straight Taper Point Needles**





Needle	Gauge	Suture Color	Length (cm)	Other information	Number per box	Code
 ST-4	19 mm	4-0	●	45	12	18507G
		4-0	●	75	12	K861G
		4-0	●	75	36	K861H





# SILK




## Braided Silk / Virgin silk










1/2 Circle TAPERCUT* Needles						
Needle	Gauge	Suture Color	Length (cm)	Other information	Number per box	Code
 V-5 17 mm	4-0	●	75		36	K952H
	3-0	●	45		12	18514G
	3-0	●	75		36	K953H
 V-7 28 mm	3-0	●	75		12	K983G
	2-0	●	75		36	K984H


3/8 Circle TAPERCUT* Needles						
Needle	Gauge	Suture Color	Length (cm)	Other information	Number per box	Code
 V-4 17 mm	4-0	●	75		36	K942H
	3-0	●	75		12	K943G

Compound Curve TAPERCUT* Needles						
Needle	Gauge	Suture Color	Length (cm)	Other information	Number per box	Code
 DA-1  VISI-BLACK 13 mm	4-0	●	45	VISI-BLACK*	12	18534G
 DA  VISI-BLACK 25 mm	4-0	●	45	VISI-BLACK*	12	18529G

1/2 Circle Reverse Cutting Needles						
Needle	Gauge	Suture Color	Length (cm)	Other information	Number per box	Code
 M-1 13 mm	5-0	●	45		36	688H
 J-1 17 mm	4-0	●	45		12	18501G
	4-0	●	45		36	734H
 X-1 22 mm	3-0	●	45		12	18502G
	3-0	●	45		36	632H
	2-0	●	45		36	737H
 CP-1 36 mm	2-0	●	75		36	443H




▼ 3/8 Circle Reverse Cutting PRIME® Needles						
Needle	Gauge	Suture Color	Length (cm)	Other information	Number per box	Code
 P-1 11 mm	6-0	●	45		36	639H
	6-0	●	45		12	786G
 P-1 P-1 11 mm	6-0	●	45		12	780G
 P-3 13 mm	5-0	●	45		36	640H
	4-0	●	45		12	18508G



▼ 3/8 Circle Reverse Cutting Needles						
Needle	Gauge	Suture Color	Length (cm)	Other information	Number per box	Code
 C-3 13 mm	5-0	●	45		36	719H
	4-0	●	45		36	736H
	4-0	●	45		12	W811
 C-3 C-3 13 mm	4-0	●	30		36	735H
	4-0	○	45		12	W816
 M-2 13 mm	4-0	●	45		36	739H
 FS-3 16 mm	6-0	●	45		12	681G
 FS-2 19 mm	5-0	●	45		36	682H
	4-0	●	45		12	EH7149G
	3-0	●	45		36	EH7345H
 FS-1 24 mm	4-0	●	45		36	629H
	3-0	●	45		12	684G
	3-0	●	45		36	684H
 FS 26 mm	2-0	●	45		12	685G
	2-0	●	45		36	685H
 FSL 30 mm	2-0	●	45		36	677H
	2-0	●	75		36	679H
	0	●	75		36	680H
 FSLX 36 mm	0	●	75		36	EH7560H



▲ 3/8 Circle Conventional Cutting PRIME® Needles						
Needle	Gauge	Suture Color	Length (cm)	Other information	Number per box	Code
 PC-3 16 mm	6-0	●	45		36	W529H
	5-0	●	45		36	W528H
	4-0	●	45		12	W604


# SILK

## Braided Silk / Virgin silk

▼ 3/8 Circle MICROPPOINT® Reverse Cutting Needles						
Needle	Gauge	Suture Color	Length (cm)	Other information	Number per box	Code
 8 mm	7-0	●	30		12	EH7464G
	6-0	●	30		12	EH7435G
 8 mm	7-0	●	45		12	768G
	6-0	●	45		12	769G
 13 mm	4-0	●	45		12	783G

▼ 1/2 Circle MICROPPOINT® Reverse Cutting Needles						
Needle	Gauge	Suture Color	Length (cm)	Other information	Number per box	Code
 8 mm	8-0	●	45		12	W870
 8 mm	8-0	●	45		12	W818

▼ 3/8 Circle MICROPPOINT® Spatula						
Needle	Gauge	Suture Color	Length (cm)	Other information	Number per box	Code
 6,5 mm	8-0	●	30		12	W1782
 6,5 mm	8-0	●	45		12	W1819

▼ 1/4 Circle Spatula Needles						
Needle	Gauge	Suture Color	Length (cm)	Other information	Number per box	Code
 8 mm	6-0	●	45		12	1780G
	5-0	●	45		12	EH7440G

Ligatures						
	Gauge	Suture Color	Length (cm)	Other information	Number per box	Code
LIGAPAK*	4-0	●	250	Reel	20	F4202
SUTUPAK*	2-0	●	140		36	EH6804H
	1	●	140		36	EH6806H
	2	●	140		36	EH6807H
	3	●	140		36	EH6808H
SUTUPAK*	3-0	●	2x70		36	EH6833H
	0	●	2x70		24	EH6835E
	1	●	2x70		24	EH6836E
	2	●	2x70		24	EH6837E
SUTUPAK*	4-0	●	3x45		36	EH6762H
	3-0	●	3x45		36	EH6763H
	2-0	●	3x45		36	EH6764H
	0	●	3x45		36	EH6765H
	2	●	3x45		36	EH6767H
	3	●	3x45		36	EH6768H
SUTUPAK*	3-0	●	5x70		36	EH6823H
	2-0	●	5x70		36	EH6824H
	0	●	5x70		36	EH6825H
SUTUPAK*	2-0	●	10x75		12	W223
	1	●	10x75		12	W225
SUTUPAK*	4-0	●	12x45		24	EH6722E
	3-0	●	12x45		24	EH6723E
	2-0	●	12x45		36	EH6724H
SUTUPAK*	4-0	●	13x60		12	W211