

PROLENE[®]

Monofilament Polypropylene

SUTURE

Description:

Monofilament	For smooth passage through tissue.
Synthetic	For minimal tissue reaction.
Non-absorbable	Provides permanent tensile strength retention in tissue.

Color:

Blue or clear.

Range:

10-0 to 1 (USP). Supplied as needled sutures, armed with, for example, VISI-BLACK[®] needles, MULTIPASS[®] needles or CC-needles. Also available with pledgets.

PROLENE[®] sutures are also available as HEMO-SEAL[®], with suture diameter to needle diameter ratio close to 1. HEMO-SEAL[®] is particularly useful in reducing needle hole leakage (for example in synthetic vascular grafts).













Indications:

PROLENE[®] sutures are intended for use in general soft tissue approximation and/or ligation, including use in cardiovascular, ophthalmic and neurosurgical procedures. Typical areas of use include:













- Cardiovascular anastomoses
- Subcuticular sutures



Full details in the Instruction for Use included in every package.














 3/8 Circle Taper Point Needles						
Needle	Gauge	Suture Color	Length (cm)	Other information	Number per box	Code
BV130-5  VISI-BLACK	6,5 mm	8-0	45	VISI-BLACK*	12	X916G
		8-0	60	VISI-BLACK*	12	KHH5663G
		8-0	75	VISI-BLACK*	24	EH7459E
		7-0	7	VISI-BLACK*	36	EH7430H
BV130-5 	6,5 mm	8-0	45	MULTIPASS*	24	EH7470E
		8-0	45	MULTIPASS*	12	W2777
		8-0	60	MULTIPASS*	36	KHH5662H
		8-0	75	MULTIPASS*	24	EH7469E
		8-0	75	Straight pack, MULTIPASS*	12	EH7469LG
		7-0	2x75	Straight pack, MULTIPASS*	12	KHL5673JG
BV130-5 	6,5 mm	8-0	13		12	W2775
BV175-6  VISI-BLACK	8 mm	8-0	60	VISI-BLACK*	36	8841H
		7-0	60	VISI-BLACK*	36	8335H
		7-0	2x60	Straight pack, VISI-BLACK*	12	M8935
BV175-6 	8 mm	8-0	60	MULTIPASS*	36	8741H
		7-0	60	MULTIPASS*	36	8735H
		7-0	2x60	Straight pack, MULTIPASS*	12	M8737G
BV175-7  VISI-BLACK		8-0	2x60	Straight pack, VISI-BLACK*	12	M8967
BV175-7 	8 mm	7-0	60		36	8766H
		7-0	2x60	Straight pack	12	M8766
BV175-8  VISI-BLACK	9,3 mm	7-0	45	VISI-BLACK*	36	EH7225H
		7-0	75	VISI-BLACK*	36	EH7222H
		7-0	2x75	Straight pack, VISI-BLACK*	12	MEH7224JG
		7-0	4x60	VISI-BLACK*	12	X967G
BV175-8 	9,3 mm	7-0	45	MULTIPASS*	36	8701H
		7-0	60	MULTIPASS*	12	W8702
		7-0	75	MULTIPASS*	36	EH7405H
		7-0	4x60	Straight pack, MULTIPASS*	12	W8702S
		7-0	4x75	Straight pack, MULTIPASS*	12	M8702G
BV175-8 	9,3 mm	7-0	75	MULTIPASS*	36	EH7402H
BV-1  VISI-BLACK	9,3 mm	7-0	60	VISI-BLACK*	12	W8304
		7-0	2x60	VISI-BLACK*, Straight pack	12	KAA5663LG
		6-0	60	VISI-BLACK*	12	W8305
		6-0	75	VISI-BLACK*	36	EH7223H





PROLENE*
Monofilament Polypropylene






3/8 Circle Taper Point Needles									
Needle	Gauge	Suture Color	Length (cm)	Other information	Number per box	Code			
 BV-1	9,3 mm	6-0	●	75		36	EH7406H		
 BV-1	9,3 mm	6-0	●	75		36	EH7403H		
 BV-1 X VISI-BLACK	9,3 mm	7-0	●	75	VISI-BLACK*, Heavier Wire Diameter	36	EH7238H		
		6-0	●	75	VISI-BLACK*, Heavier Wire Diameter	36	EH7240H		
 BV-1 X	9,3 mm	6-0	●	75	Heavier Wire Diameter	36	EH7245H		
 BV	11 mm	6-0	●	75		36	8776H		
 C-1 VISI-BLACK	13 mm	6-0	●	45	VISI-BLACK*	36	X577H		
		6-0	●	60	VISI-BLACK*	12	W8307		
		6-0	●	75	VISI-BLACK*	36	EH7228H		
		6-0	●	90	VISI-BLACK*	36	EH7230H		
		6-0	●	2x75	Straight pack, VISI-BLACK*	12	MEH7232JG		
		5-0	●	75	VISI-BLACK*	36	EH7229H		
		5-0	●	90	VISI-BLACK*	36	EH7231H		
		5-0	●	2x75	Straight pack, VISI-BLACK*	12	MEH7233JG		
 C-1	13 mm	7-0	●	60		12	W8725		
		7-0	●	75		36	8708H		
		6-0	●	60		12	W8718		
		6-0	●	75		36	8708H		
		6-0	●	2x90	Straight pack	12	KBL5690JG		
		5-0	●	45		36	8717H		
		5-0	●	60		36	KBB5661H		
		5-0	●	70	Pledget	12	KBB5671VG		
		5-0	●	75		36	EH7477H		
		5-0	●	90		36	8720H		
		5-0	●	2x90	Straight pack	12	KBL5691JG		
		4-0	●	75		36	EH7478H		
		 C-1 VISI-BLACK	13 mm	6-0	●	75	VISI-BLACK*	36	EH7226H
				5-0	●	75	VISI-BLACK*	36	EH7227H
 C-1	13 mm	6-0	●	75		36	8889H		
		5-0	●	75		36	8890H		
		4-0	●	75		36	EH7476H		
 C-1 X VISI-BLACK	13 mm	5-0	●	75	VISI-BLACK*, Heavier Wire Diameter	36	EH7257H		
 C-1 X	13 mm	5-0	●	75	Heavier Wire Diameter	36	EH7255H		
 C-1 X	13 mm	5-0	●	75	Heavier Wire Diameter	36	EH7254H		

3/8 Circle Taper Point Needles						
Needle	Gauge	Suture Color	Length (cm)	Other information	Number per box	Code
 BB	17 mm	5-0	●	90		36 8580H
		4-0	●	90		36 8581H
 BB	17 mm	5-0	●	75		36 8880H
		4-0	●	75		36 8881H



1/2 Circle Taper Point Needles						
Needle	Gauge	Suture Color	Length (cm)	Other information	Number per box	Code
 TF-6	6,5 mm	7-0	●	45		24 EH7813E
		7-0	●	75		36 PQQ5671H
		6-0	●	45		24 EH7814E
		6-0	●	60		24 PQQ5660E
 TF-1	10 mm	6-0	●	75		36 EH7835H
		5-0	●	60		36 PUU5662SH
 RB-2	13 mm	6-0	●	45		36 8714H
		6-0	●	75		36 8711H
		5-0	●	75		36 8710H
		5-0	●	4x75	Straight pack	12
 TF	13 mm	6-0	●	75		36 EH7400H
		5-0	●	60		36 8205H
		5-0	●	75		36 EH7401H
		4-0	●	90		36 EH7411H
 RB-1 VISI-BLACK	17 mm	5-0	●	90	VISI-BLACK*	36 VB50H
		4-0	●	90	VISI-BLACK*	36 8357H
 RB-1	17 mm	5-0	●	90		36 8556H
		4-0	●	90		36 8557H
		3-0	●	90		36 8558H
 RB-1	17 mm	5-0	●	75		36 8870H
		4-0	●	75		36 8871H
		3-0	●	75		36 8872H
 JB-1 VISI-BLACK	22 mm	4-0	●	90	Pledgets, VISI-BLACK*	36 EH8065H
 SH-1	22 mm	4-0	●	90		36 EH7585H
		3-0	●	90		36 EH7586H
 SH-1	22 mm	4-0	●	75		36 EH7581H
		3-0	●	75		36 EH7582H
 JB VISI-BLACK	26 mm	4-0	●	90	VISI-BLACK*	36 EH8087H
		4-0	●	90	Pledgets, VISI-BLACK*	36 EH8062H
		4-0	●	120	VISI-BLACK*	36 X1166H
		3-0	●	90	VISI-BLACK*	36 EH8086H

PROLENE*
Monofilament Polypropylene


● 1/2 Circle Taper Point Needles						
Needle	Gauge	Suture Color	Length (cm)	Other information	Number per box	Code
 SH 26 mm	4-0	●	90		36	8521H
	4-0	●	90	Pledgets, Straight pack	12	EH7190LG
	3-0	●	90		36	8522H
	3-0	●	120		36	EH7584H
	2-0	●	90		36	8523H
 SH 26 mm	4-0	●	75		36	8831H
	3-0	●	75		36	8832H
	2-0	●	75		36	8833H
	0	●	75		36	8834H
 JMH-1 31 mm VISI-BLACK	3-0	●	120		36	X1295H
 MH 36 mm	3-0	●	90		36	EH7499H
	3-0	●	120		36	PGG5613H
	2-0	●	120		36	PGG5614H
	0	●	90		36	PGG5695H

● 1/2 Circle Taper Point Needles (Heavier Wire Diameter)						
Needle	Gauge	Suture Color	Length (cm)	Other information	Number per box	Code
 CT-2 26 mm	2-0	●	75		36	8411H
	0	●	75		36	8412H
	1	●	90		36	PK5696H
 CT-1 36 mm	3-0	●	75		36	8422H
	2-0	●	75		36	8423H
	0	●	75		36	8424H
	1	●	75		36	8425H
 CT 40 mm	0	●	75		36	8434H
	1	●	100		12	W742
 CTX 48 mm	0	●	90		36	PN5695H
	1	●	150	Loop	12	8527G
 CTX 48 mm	0	●	90		24	PNN5695E




Compound Curve Taper Point Needles

Needle	Gauge	Suture Color	Length (cm)	Other information	Number per box	Code
 ABV-1	9,3 mm	7-0	75		36	EH8000H
 AC-1	13 mm	6-0	75		36	EH8010H
		5-0	75		36	EH8014H







5/8 Circle Taper Point Needles

Needle	Gauge	Suture Color	Length (cm)	Other information	Number per box	Code
 UR-6	26 mm	2-0	75		36	VC5673H




Straight Taper Point Needles








Needle	Gauge	Suture Color	Length (cm)	Other information	Number per box	Code
 ST-4	19 mm	4-0	35		36	EH7508H
 ST-1	45 mm	3-0	35		36	EH7921H
 ST-79	70 mm	2-0	75		12	W8400






3/8 Circle CC-Needles





Needle	Gauge	Suture Color	Length (cm)	Other information	Number per box	Code
 CC VISI-BLACK	8 mm	8-0	60	VISI-BLACK*	12	W8101
 CC	9,3 mm	8-0	60		24	EH7972E
		8-0	60	MULTIPASS*	12	W8703
		7-0	45		36	8705H
		7-0	75	MULTIPASS*	36	8704H
		6-0	60	MULTIPASS*	12	W8807
		5-0	60	HEMOSEAL*	12	W8862
 CC	9,3 mm	7-0	75		36	EH7970H
 CCX VISI-BLACK	9,3 mm	7-0	75	VISI-BLACK*, Heavier Wire Diameter	36	EH7237H
		6-0	75	VISI-BLACK*, Heavier Wire Diameter	36	EH7234H
 CCX	9,3 mm	7-0	75	Heavier Wire Diameter	36	EH7236H
		6-0	75	Heavier Wire Diameter	36	EH7235H
 CC-11	11 mm	6-0	60		12	W8802







PROLENE*
Monofilament Polypropylene

3/8 Circle CC-Needles						
Needle	Gauge	Suture Color	Length (cm)	Other information	Number per box	Code
 CC-1 13 mm	8-0	●	45		36	8719H
	6-0	●	75	MULTIPASS*	36	8707H
	6-0	●	2x75	MULTIPASS*, Straight pack	12	X535G
	5-0	●	75		36	EH7977H
	5-0	●	90		36	8721H
 CC-1 13 mm	6-0	●	75		36	EH7975H
	5-0	●	75		36	EH7976H
 CC-1 X 13 mm	5-0	●	75	Heavier Wire Diameter	36	EH7251H



3/8 Circle TAPERCUT* Needle						
Needle	Gauge	Suture Color	Length (cm)	Other information	Number per box	Code
 MW-16, MW-14 5,5 mm	8-0	●	60		24	EH7828E
 V-18, V-18 13 mm	6-0	●	75		36	EH7474H
	5-0	●	75		36	EH7475H
 V-18 13 mm	6-0	●	75		36	EH7471H
	5-0	●	75		36	EH7472H
 V-18 X, V-18 X 13 mm	5-0	●	75	Heavier Wire Diameter	36	EH7253H
 V-4 17 mm	5-0	●	75		36	F2829
	4-0	●	75		36	F2828
	4-0	●	90		36	8925H
 V-4 17 mm	5-0	●	75		36	8941H
	4-0	●	75		36	8942H
	3-0	●	75		36	8943H
 V-26 26 mm	3-0	●	90		36	EH7788H



1/2 Circle CC-Needles						
Needle	Gauge	Suture Color	Length (cm)	Other information	Number per box	Code
 CC-6 11 mm VISI-BLACK	7-0	●	60	VISI-BLACK*	12	X1004G
 CC-4 13 mm VISI-BLACK	6-0	●	75	VISI-BLACK*	36	EH7242H
	5-0	●	75	VISI-BLACK*	36	EH7241H
	5-0	●	2x75	VISI-BLACK*, Straight pack	12	MEH7241JG
 CC-4 13 mm	6-0	●	60	MULTIPASS*	12	W8815
	5-0	●	75		12	W8816
	5-0	●	60	HEMO-SEAL*	12	W8664
	4-0	●	75	HEMO-SEAL*	12	W8665
 CC-16 16 mm	3-0	●	90	HEMO-SEAL*	12	W8667
 CC-30 31 mm	3-0	●	90		12	W8849


 Compound Curve CC-Needles						
Needle	Gauge	Suture Color	Length (cm)	Other information	Number per box	Code
 VISI-BLACK	9,3 mm	7-0	60	VISI-BLACK*	36	EH8020H
		7-0	75	VISI-BLACK*	24	EH8021E
		6-0	75	VISI-BLACK*	36	EH8025H
 VISI-BLACK	13 mm	6-0	75	VISI-BLACK*	36	EH8030H
		6-0	90	VISI-BLACK*	36	EH8031H
		6-0	2x75	Straight pack	12	MEH8034JG
		5-0	75	VISI-BLACK*	36	EH8036H
		5-0	90	VISI-BLACK*	36	EH8037H
 VISI-BLACK	13 mm	6-0	75	VISI-BLACK*	36	EH8032H





 1/2 Circle TAPERCUT* Needles						
Needle	Gauge	Suture Color	Length (cm)	Other information	Number per box	Code
	17 mm	5-0	90		36	8934H
		4-0	90		36	8935H
		4-0	90	Straight pack, Pledgets	12	8954G
		4-0	2x90	Straight pack	12	OAL5693JG
		3-0	90		36	8936H
		2-0	90		36	8937H
	17 mm	5-0	75		36	8951H
		4-0	75		36	8952H
		3-0	75		36	8953H
	22 mm	4-0	90		36	EH7774H
		3-0	90		36	EH7775H
	22 mm	2-0	75		36	EH7772H
	26 mm	4-0	90		36	8975H
		4-0	120		36	OBB5613H
		3-0	90		36	8976H
		3-0	90	Pledgets	36	OBB5694SH
		3-0	2x90	Straight pack	12	OBL5694JG
		3-0	120		36	EH7778H
		2-0	90		36	8977H
		2-0	4x90	Straight pack, Pledgets	12	OBL5695LG
		2-0	4x90	Straight pack	12	OBL5695LG
		2-0	120		36	EH7779H
		0	90		36	8978H
		0	120		36	OBB5616H






PROLENE*
Monofilament Polypropylene





1/2 Circle TAPERCUT* Needles						
Needle	Gauge	Suture Color	Length (cm)	Other information	Number per box	Code
 V-7 26 mm	4-0	●	75		36	8962H
	3-0	●	75		36	8963H
	2-0	●	75		36	8964H
 V-34 36 mm	2-0	●	75		36	8443H
	0	●	75		36	8444H

Compound Curve TAPERCUT* Needles						
Needle	Gauge	Suture Color	Length (cm)	Other information	Number per box	Code
 DA-3 VISI-BLACK 9,3 mm	7-0	●	45	VISI-BLACK*	12	18541G
 DA-2 VISI-BLACK 13 mm	6-0	●	45	VISI-BLACK*	12	18543G





Straight TAPERCUT* Needles						
Needle	Gauge	Suture Color	Length (cm)	Other information	Number per box	Code
 V-47 38 mm	4-0	●	35		36	8592H




3/8 Circle Reverse Cutting PRIME* Needles						
Needle	Gauge	Suture Color	Length (cm)	Other information	Number per box	Code
 P-6 8 mm	7-0	●	45		36	8648H
	6-0	●	45		24	W8000T
 P-1 11 mm	7-0	●	75	MULTIPASS*	36	MPP2854H
	6-0	●	45		36	8697H
	6-0	●	45	MULTIPASS*	36	MPP8697H
 P-3 13 mm	6-0	●	45		36	8695H
	6-0	●	45	MULTIPASS*	36	MPP8695H
	5-0	●	45		36	8698H
	5-0	●	45	MULTIPASS*	36	MPP8698H
	5-0	○	75	MULTIPASS*	36	MPP2858H
	4-0	●	45		36	8699H
	4-0	●	45	MULTIPASS*	36	MPP8699H
	4-0	○	45		36	JA5742H
	4-0	○	45	MULTIPASS*	36	MPP8604H
	 PS-3 16 mm	6-0	●	75	MULTIPASS*	36
5-0		●	45		36	8681H
5-0		●	45	MULTIPASS*	36	MPP8681H
4-0		●	45		36	EH7764H
4-0		●	45	MULTIPASS*	36	MPP7764H


3/8 Circle Reverse Cutting PRIME* Needles						
Needle	Gauge	Suture Color	Length (cm)	Other information	Number per box	Code
 PS-2 19 mm	5-0	●	45		24	W8014T
	5-0	●	45	MULTIPASS*	36	MPP8686H
	4-0	●	45		24	W8015T
	4-0	●	45	MULTIPASS*	36	MPP8682H
	4-0	○	75	MULTIPASS*	36	MPP2861H
	3-0	●	60	MULTIPASS*	36	MPP2835H
	3-0	●	75		36	F2835
	3-0	○	75	MULTIPASS*	36	MPP2862H
 PS-2 S 19 mm	4-0	●	45	Slimmed	36	8682H
 PS-1 24 mm	4-0	●	45		36	EH7758H
	3-0	●	45	MULTIPASS*	36	MPP7757H
	3-0	●	60		12	EH7666G
	3-0	●	75		36	EH7568H
 PS 26 mm	2-0	●	45	MULTIPASS*	36	MPP8026H
 PSL 30 mm	2-0	●	75		36	EH7121H


3/8 Circle Reverse Cutting Needles						
Needle	Gauge	Suture Color	Length (cm)	Other information	Number per box	Code
 M-3 VISI-BLACK 16 mm	6-0	●	45	VISI-BLACK*	12	8750G
 FS-3 16 mm	6-0	●	45		36	8660H
	5-0	●	45		36	8614H
 FS-2 19 mm	5-0	●	45		36	8661H
	4-0	●	45		36	EH7151H
	4-0	●	75		36	EH7150H
	3-0	●	45		36	8665H
	3-0	●	75		36	EH7693H
 FS-1 24 mm	4-0	●	45		36	8629H
	3-0	●	45		36	8694H
	3-0	●	60		12	EH7963G
	3-0	●	75		36	EH7694H



PROLENE*
Monofilament Polypropylene








▼ 3/8 Circle Reverse Cutting Needles						
Needle	Gauge	Suture Color	Length (cm)	Other information	Number per box	Code
 FS 26 mm	2-0	●	45		36	8685H
 FSL 30 mm	4-0	●	75		36	8670H
	3-0	●	75		36	8675H
	2-0	●	75		36	EH7697H
	0	●	45		36	EH7920H
 FSLX 36 mm	3-0	●	75		36	8649H
	2-0	●	75		36	8689H
	0	●	75		36	EH7561H
 LR 76 mm	2	●	3x50		24	D9802


▲ 3/8 Circle Conventional Cutting PRIME* Needles						
Needle	Gauge	Suture Color	Length (cm)	Other information	Number per box	Code
 PC-3 16 mm	6-0	●	45		24	W8873T
	6-0	●	45	MULTIPASS*	36	MPP8873H
	5-0	●	45		12	W527
	5-0	●	45		24	W8874T
	5-0	●	45	MULTIPASS*	36	MPP8874H
	4-0	●	45		12	W525
	4-0	●	45		24	W8875T
	4-0	●	45	MULTIPASS*	36	MPP8875H
 PC-5 19 mm	5-0	●	45	MULTIPASS*	36	MPP8876H
	4-0	●	45		24	W8877T
	4-0	●	45	MULTIPASS*	36	MPP8877H
 PC-25 26 mm	4-0	●	45		12	W539
	3-0	●	45		12	W538




▼ 1/2 Circle Reverse Cutting PRIME* Needles						
Needle	Gauge	Suture Color	Length (cm)	Other information	Number per box	Code
 PS-4 16 mm	4-0	●	45		36	EH7759H

▼ 1/2 Circle Reverse Cutting Needles						
Needle	Gauge	Suture Color	Length (cm)	Other information	Number per box	Code
 J-1 17 mm	3-0	●	45		36	EH7545H


▼ Straight Cutting Needles							
Needle	Gauge	Suture Color	Length (cm)	Other information	Number per box	Code	
 SC-1	13 mm	4-0	●	45		36	EH7730H
 KS	60 mm	3-0	●	75		36	8622H
		2-0	●	75		36	8623H


● 3/8 Circle Taper Point Needles Special							
Needle	Gauge	Suture Color	Length (cm)	Other information	Number per box	Code	
 BV75-3	3,8 mm	10-0	●	13		12	W2790
 BV130-4	4,7 mm	9-0	●	13		12	W2780
		8-0	●	13		12	W2779
 BV-6	4,8 mm	10-0	●	10		12	EH7595G
 BV-4	4,9 mm	10-0	●	15		12	EH7596G
 BV100-4	5,1 mm	10-0	●	13		12	W2794
		9-0	●	13		12	W2783
 BV100-4	5,1 mm	9-0	●	45		12	F1891
 BV-2	6,5 mm	9-0	●	15		12	EH7597G



▼ 1/4 Circle TAPERCUT* Needles							
Needle	Gauge	Suture Color	Length (cm)	Other information	Number per box	Code	
 CIF-4	13 mm	10-0	●	20		12	788G





▼ 3/8 Circle MICROPOINT* Spatula							
Needle	Gauge	Suture Color	Length (cm)	Other information	Number per box	Code	
 GS-12	5,5 mm	10-0	●	30		12	EH7473G
 GS-12 GS-12	5,5 mm	10-0	●	30		12	EH7598G
 CS-35-C CS-35-C	6 mm	10-0	●	30		12	EH7994G

PROLENE*
Monofilament Polypropylene

1/4 Circle CS ULTIMA* Spatula							
Needle	Gauge	Suture Color	Length (cm)	Other information	Number per box	Code	
 CTC-6 L 13 mm	10-0	●	10		12	EH7897G	

Ski MICROPOINT* Spatula							
Needle	Gauge	Suture Color	Length (cm)	Other information	Number per box	Code	
 CSS-160 CSS-160 16 mm	10-0	●	30		12	EH7997G	

Straight MICROPOINT* Spatula							
Needle	Gauge	Suture Color	Length (cm)	Other information	Number per box	Code	
 STC-6 STC-6 16 mm	10-0	●	20		12	1713G	
 STC-6 16 mm	10-0	●	15	Loop	12	EH7900G	

Special Combinations							
Needle	Gauge	Suture Color	Length (cm)	Other information	Number per box	Code	
 STC-6 CS-35-C 16 mm + 6 mm	10-0	●	20		12	ET5631SG	
 GS-9 STC-6 6,5 mm + 16 mm	9-0	●	20		12	1753G	
 ST-1 ST-1 FS 45 mm + 45 mm + 26 mm	3-0	●	35	Loop with polypropylene button. Pull-out tendon suture for Bunnell technique.	12	EH7922G	
 V-47 V-47 FS-1 38 mm + 38 mm + 26 mm	4-0	●	35	Loop with polypropylene button. Pull-out tendon suture for Bunnell technique.	12	EH7668G	