

# PRODUCT DATA SHEET

INFUSOMAT SP. LINE, STAND., PVC, LL, 250CM (8700435SP)



Attribute Name	Attribute Value
Brand name	Infusomat®
Product name (HC)	Infusomat® Space Line
Product type	Administration set
Application	IV-Standard
Type of administration set	Pressure
Container connector	IV standard spike
Length	250.00 cm
Length from pump outlet	145 cm
Inner-ø tube	3.00 mm
Outer-ø tube	4.10 mm
Type of filter	Fluid filter
Position of filter	Drip chamber
Pore size (Fluid filter)	15 µm
drops/ml	20
Air ventilation	vented

Attribute Name	Attribute Value
Flow control device	Roller clamp
Injection access	-
with AirStop	No
with PrimeStop	No
Number of back check valves	-
Number of SafSite®	-
Number of needle-free valves	-
Connection	Luer-Lock/male
Light-protected	No
Color of light protection	-
Pump compatibility (Type)	Infusomat® fm/fmS, Infusomat® Space
Length of pump segment	100.00 mm
max. pressure	2.0 bar
Material	PVC
Catalogue attribute	IV-Standard
Type of packaging	Ward pack
Product (unpack.) contains PVC	Yes
Product (unpack.) cont. Latex	No
Product (unpack.) cont. DEHP	No
Packaging contains PVC	No
Packaging contains Latex	No
Packaging contains DEHP	No
Product contains PVC	Yes
Product contains Latex	No
Product contains DEHP	No
Product contains Nickel	No
Product contains BPA	No
Disposable / Reusable	disposable
Sterility	sterile
Type of Sterilization	E0
Type of Sterilizat. (allowed)	E0

Attribute Name	Attribute Value
UMDNS-Version	1.0
UMDNS	16579-Infusion Pump Adminis...
Priming volume	17.400 ml
Length (Inch)	98
Type of Injection Access	-
eCl@ss Code (Rel. 11)	34220605
eCl@ss Description (Rel. 11)	Inf. pump tub. (electr. power)