

<b>Product Name</b>	Paediatric respi-check high concentration mask		
<b>Reference numbers</b>	1292000		
<b>Product Type</b>	Respi-check mask		
<b>Manufacturer</b>	Intersurgical		
<b>Materials</b>	PVC PP	Aluminium Polystyrene ball	Elastic
<b>Packaging</b>	Individually packaged		
<b>Unit Packaging</b>	20		
<b>Box quantity</b>	640		
<b>Pallet quantity</b>			
<b>Instructions for Use</b>	Instructions for use with product	YES	
	Instructions for use with each box	YES	
<b>Certification</b>	CE0120	Class 2A	
<b>Applicable Standards</b>	BS EN ISO 13485:2012 Medical devices — Quality management systems — Requirements for regulatory purposes		
	ISO 9001:2008 Quality Management Systems - Requirements		
	ISO 15223-1:2012 Medical Devices - Symbols to be Used with Medical Device Labels, Labelling, & Information to be Supplied		
	BS EN 13544-2:2002 Respiratory therapy equipment - Tubing and connectors		
<b>Sterile</b>	No		
<b>Shelf life</b>	5 years		
<b>Disposal Methods</b>	As hospital guidelines		
<b>Re-use</b>	Single use only		



**Product characteristics**

**Mask information**

Height (mm)	97
Width (mm)	73
Approx dead space (ml)	90

**Tubing information**

Weight (g)	37
Length (m)	2.1
Outside diameter (mm)	5.6

**Tubing has a star shaped lumen**

**Bag information**

Height (mm)	229
Width (mm)	156
Capacity (l)	1

**Total weight of mask + bag = 36.2**

The Paediatric Respi-Check High Concentration mask allows continuous monitoring of the respiratory rate whilst delivering oxygen to the patient. It can be connected to a Respi-Check monitor which accurately records the respiratory rate.

