






Dynacast Prelude - Gipsskena		Best.nr.	Storlek	Antal/fp	Antal/trp
Glasfiber-gipsskena med polster. Kardborrband kit.			Vit		
Allt-i-ett-system.		71418-00	2,5 cm x 4,6 m 5-lagers	1	2
Latexfri.		71420-00	5 cm x 4,6 m 5-lagers	1	2
		71421-00	7,5 cm x 4,6 m 7-lagers	1	2
		71415-00	10 cm x 4,6 m 7-lagers	1	2
		71416-00	12,5 cm x 4,6 m 7-lagers	1	2
		71417-00	15 cm x 4,6 m 7-lagers	1	2
		71419-00	20 cm x 2,3 m 7-lagers	1	2
			Röd		
		71415-18	1,5cm x 42,5 1st 5cm x 49cm 2st	10	10
			Blå		
		71415-19	1,5cm x 42,5 1st 5cm x 49cm 2st Allt-ett-system med polster	10	10