



## Y'-Extension set 2-way

### MFx2291

- Set Length 209 cm
- Priming Volume 2.3 ml
- Bore 1 mm
- 2 Way Y Connector
- 1 Pinch Clamp(s)
- 1 Back-Check Valve(s)
- 1 Female Luer(s)
- 1 Anti-Syphon Valve
- 1 Male Luer(s)
- Rotating Male Luer(s)
- Natural rubber latex is not part of the material formulation
- DEHP is not part of the material formulation
- Quantity 50

CE