

MF2292E



# Total Intravenous Administration (TIVA) Extension Set



**SKU/REF** MF2292E

**SKU/REF Name** TIVA extension set 3 way

## Description

2 Anti-Siphon Valves - 1 Back Check Valve - BD SmartSite™ needle-free connector - length 210 cm

### Resources

Resources (1)

 [Infusion specialty disposables e-catalogue](#)

## Specifications

<b>Country of Manufacture</b>	Malta
<b>Classification</b>	IIa
<b>Quantity per Box</b>	150 (50x3)
<b>Tube Length (cm)</b>	206
<b>Total Length (cm)</b>	210
<b>Priming Volume (mL)</b>	3.9
<b>Tube Inner Diameter (mm)</b>	1 (second tube 3)
<b>Pump Use</b>	Yes, maximum 2.5 bar
<b>Shelf Life</b>	5 Years
<b>Sterilization Method</b>	Ethylene Oxide
<b>Polyvinyl chloride (PVC)</b>	Yes
<b>DEHP</b>	DEHP-Free
<b>Ways</b>	3-Way
<b>Back check valve</b>	Yes
<b>Anti-syphon valve</b>	Yes
<b>Luer Lock</b>	Rotating luer lock
<b>BD SmartSite™</b>	Yes
<b>Application Area</b>	Anaesthesia

