





# PERFORMANCE UNDER PRESSURE



## PRESSURE INFUSION BAG

**Dual pressure valve:** Standard pressure all the time, over-pressurize when needed with the click of a button.



# Merit Medical Pressure Infusion Bag

- Two pressure settings
- Available in 500ml, 1000ml, and 3000ml sizes
- Loads quickly and easily
- Non-Sterile



## STANDARD PRESSURE

Valve is factory pre-set for pressurization up to 300mmHg



## OVER PRESSURE

Depress blue button to set valve for over-pressure mode allowing 300-450mmHg for rapid infusion



## POLE HANGER

Easy and stable placement onto IV pole



## FLUID BAG HOOK

Bag attaches easily to rigid hook which secures bag in the event of pressure loss



## ONE-HANDED PREP

Stopcock placed close to bulb for easy inflation, closure once fully pressurized, and deflation

## ORDERING INFORMATION

Product Number	Description	Quantity
PIB500	500ml Infusor Bag with Dual Pressure Valve	12 per box, 48 per case
PIB1000	1000ml Infusor Bag with Dual Pressure Valve	12 per box, 48 per case
PIB3000	3000ml Infusor Bag with Dual Pressure Valve	6 per box, 24 per case

Before using, refer to Instructions for Use for indications, contraindications, warnings, precautions and directions for use.



Understand. Innovate. Deliver.™

Merit Medical Systems, Inc.  
1600 West Merit Parkway  
South Jordan, Utah 84095  
1.801.253.1600  
1.800.35.MERIT

Merit Medical Europe, Middle East, & Africa (EMEA)  
Amerikalaan 42, 6199 AE Maastricht-Airport  
The Netherlands  
+31 43 358 82 22

Merit Medical Ireland Ltd.  
Parkmore Business Park West  
Galway, Ireland  
+353 (0) 91 703 733

merit.com

Australia  
+61 (0) 3 8373 4943

Austria  
0800 295 374

Belgium  
0800 72 906 (Dutch)  
0800 73 172 (Français)

Brazil  
+55 11 4561 - 0062

Canada  
1.800.364.4370

China  
+86-10 8561 0788

Denmark  
80 88 00 24

Finland  
0800 770 586

France  
0800 91 60 30

Germany  
0800 182 0871

Hong Kong  
+852-2207 6888

India  
+91-80-41223376

Ireland (Republic)  
1800 553 163

Italy  
800 897 005

Japan  
800 897 005

Luxembourg  
8002 25 22

Netherlands  
0800 022 81 84

Norway  
800 11629

Portugal  
308 801 034

Russia  
+7 495 221 89 02

Singapore  
+65 67536311

South Korea  
+8210 2205 0707

Spain  
+34 911238406

Sweden  
020 792 445

Switzerland  
(Deutsch)  
+41 22 518 02 30

(Français)  
+41 22 518 02 52

(Italiano)  
+41 22 518 00 35

UK  
0800 973 115