



KATETERISERINGSSET

ART.NR 605030N

	Beskrivning	Avd/Transport Box
		10/40
1	Omslagsduk 34 g, 50 x 80 cm	
1	Rondskål 700 ml	
1	Skål 160 ml	
2	Kompresser Nonwoven 10 x 10 cm	
5	Tork Nonwoven Ø 50 mm	
1	Peang blå 12 cm	
1	Vikt Duk 34 g 50 x 80 cm	

Material:

Fields: 100 % Cellulose tissue + PE Polyethylene

Kidney dish: Polystyrene PS

Bowl: Polypropylene PP

Balls/Swabs: 70 % Rayon 30 % Polyester Peroxide bleached

Peang: Polypropylene PP

Packaging: Peel-pack, Paper/PE

Sterilization: Shelf-life 5 years from production date

Manufacturer: Medioplast AB

MDD Classification: Is

To be stored in room temperature

Transport Box = minimum sales unit

Destruction: Combustion, recyclable carton, recyclable well-cardboard



STERILE | EO

