

Transferpipetter



- Tillverkad i Polyetylen (PE)
- För säker överföring

Artikelnummer	Benämning	Längd mm	Dropp stl. μ l	Inner/ytter frp
86.1170	Transferpipett 1,0 mL, graderad till 1 mL (1:0.10)	115	45-55	500/2000
86.1171	Transferpipett 3,5 mL, graderad till 3 mL (3:0.50)	155	35-55	500/2000
86.1172	Transferpipett 3,5 mL, graderad till 1 mL (1:0.25)	156	30-45	500/2000
86.1173	Transferpipett 3,5 mL	155	30-45	500/2000
86.1174	Transferpipett 6,0 mL	152	40-55	400/1600
86.1175	Transferpipett 6,0 mL	146	20-30	400/1600
86.1176	Transferpipett 2,0 mL	154	35-45	750/3000
86.1179	Transferpipett 1,0 mL	87	30-45	750/3000
86.1180	Transferpipett 1,0 mL	104	17-22	750/3000