


Product Data Sheet

Kapillärrör End-to- End

REF 165713



Tryckt: 8.1.2015

REF	Storlek	
165713	25 µL	
CE IVD ISO 7750		
Material:	Glas HGB 3	
Egenskaper:	End-to-End Mikro pipetter 25µl ofärgad Physical Dimensions Internal Diameter: 0,80 mm Tolerance: ± 0,020 mm External Diameter: 1,05 mm Tolerance: ± 0,05 mm Length: 50 mm Tolerance: ± 1 mm Nominal Volume: 25 µl Tolerance: ± Acc: ≤ 0,50 % CV: ≤ 0,60 % Material: Glass, HGB 3	
Steril produkt:	Nej	
Hållbarhetstid:	4 år	
MBA/st:	100 st	
Avdelningsförpackning/st :	100 st	
Transportlåda/st:	1000 st	