

Product Data Sheet

Kapillärrör End-to- End

REF 165613



Tryckt: 3.11.2014

REF	Storlek																	
165613	20 µL																	
CE IVD ISO 7750																		
Material:	Glas HGB 3																	
Egenskaper:	<p>End-to-End Mikro pipetter 20µl svart färgkod</p> <p>Physical Dimensions</p> <table><tr><td>Internal Diameter:</td><td>0,80 mm</td><td>Tolerance: ±</td><td>0,020 mm</td></tr><tr><td>External Diameter:</td><td>1,05 mm</td><td>Tolerance: ±</td><td>0,05 mm</td></tr><tr><td>Length:</td><td>40 mm</td><td>Tolerance: ±</td><td>1 mm</td></tr><tr><td>Nominal Volume:</td><td>20 µl</td><td>Tolerance: ±</td><td>Acc: ≤ 0,50 % CV: ≤ 0,60 %</td></tr></table> <p>Material: Glass, HGB 3</p>		Internal Diameter:	0,80 mm	Tolerance: ±	0,020 mm	External Diameter:	1,05 mm	Tolerance: ±	0,05 mm	Length:	40 mm	Tolerance: ±	1 mm	Nominal Volume:	20 µl	Tolerance: ±	Acc: ≤ 0,50 % CV: ≤ 0,60 %
Internal Diameter:	0,80 mm	Tolerance: ±	0,020 mm															
External Diameter:	1,05 mm	Tolerance: ±	0,05 mm															
Length:	40 mm	Tolerance: ±	1 mm															
Nominal Volume:	20 µl	Tolerance: ±	Acc: ≤ 0,50 % CV: ≤ 0,60 %															
Steril produkt:	Nej																	
Hållbarhetstid:	4 år																	
MBA/st:	100 st																	
Avdelningsförpackning/st :	100 st																	
Transportlåda/st:	1000 st																	