




Medline International France SAS

5 rue Charles Lindbergh  
44110 Châteaubriant

[www.medline.com/uk](http://www.medline.com/uk)



## Technical sheet

# Smart IQ Breathable Impervious Surgical Gown

### Reference(s)

LW89015CEA, LW89045CEA, LW39079CEA

### Description

A sterile wrap 90x90 cm contains:  
- 1 breathable impervious surgical gown  
- 2 hand towels, 39x30 cm unfolded

Reference	Length (cm)	Size	Unit/case
LW89015CEA	120	L	24
LW89045CEA	130	XL	20
LW39079CEA	150	2XLL	20

### Main materials

Breathable impervious material (polypropylene + polyethylene)

### Packaging

Vented sterilization pouch  
2 polyliners contain all the units in a carton box

### Manufacturer

Medline International France SAS  
EN ISO 13485

### Regulation

In compliance with Medical Devices Directive 93/42/EEC  
MDD Class: I sterile  
HP: High performance according to EN 13795 standard



Latex free



0123



STERILE EO



© 2011 Medline Industries Inc. or one of its subsidiaries.  
All rights reserved.

LIT-LW95030CE-EN02

01.05.SUB.04