

# BD SmartSite™ needle-free connector

**SKU/REF** 2000E7D  
**SKU/REF Name** BD SmartSite™ needle-free connector



## Specifications

Length (cm)	2.5
Priming Volume (mL)	0.1
Maximum Pressure Rating	325 PSI@10 mL per second
Gravity Flow Rate	>6L/hr
Needle-Free Connector	BD SmartSite™
Change Interval	7days or 200 activations
Type of Luer	Male Spin Lock
DEHP	DEHP-Free
Sterilization Method	Ethylene Oxide
Country of Manufacture	Mexico
Classification	Ila
Quantity per Box	50
Quantity per Case	100
Shelf Life	3 Years

CE 2797

<https://www.bd.com/en-eu/bd-smartsite-needle-free-connector/2000e7d?return=true>

Becton Dickinson Sverige AB, Marieviksgatan 25,11743 Stockholm  
bdsweden@bd.com

**bd.com**

BD, BD-logo, NeutroX, NeutraClear är varumärken som tillhör Becton, Dickinson and Company eller dess dotterbolag  
© 2022 BD. Alla rättigheter förbehållna. BD-59783 March 2022

