



Alaris™ SE Primary infusion set 72303

- Set Length 254 cm
- Priming Volume 24 ml
- Bore 2.7 mm
- 1 Spiked Drip Chamber(s) (Vented/Unvented)
- Drops per ml 20
- Hydrophobic Air Inlet Filter (with Cap)
- 15 µm Drip Chamber Filter
- AccuSlide™ Flow Regulator
- 1 Male Luer(s)
- Polyvinyl Chloride (PVC)
- Natural rubber latex is not part of the material formulation
- Quantity 50

