


60613 HIP SET



Included Article Numbers	- 4 Cellulose Towel 18x25cm 70504 1 Split Sheet 200x260cm, adh. split 20x102cm, patched 76410 1 Adhesive OP-sheet 300x175 cm, patched 95190 1 Reinforced Table Cover 150x190cm, abs. 75x190cm 90190 1 Table Cover 150x190 cm, abs. 75x190 cm 911500 2 OP-Tape 9x49cm
Other information	Removable label on sterility barrier
Country Of Origin	Belgium
Dispenser Box Quantity	5
Transport Box Quantity	10
Standard	EN 13795-3 High Performance EN 13795-1,2,3 ISO 11607-1 ISO 10993 ISO 14001
Label Of Standard	EN 980 CEE 93/42 ISO 15223
Instruction Intended Use	Surgical drapes, when sterilised, are intended to minimize the spread of microorganisms in order to reduce the risk for post operative wound infection.
Sterilization Method	Irradiation
MDD Classification	Class I Sterile
CEMark Certificate	01966