

Blue Line Ultra® Tracheostomy Tube Kit

- * Röntgentät termosensitivt PVC material
- * 105° cirkelformad kurvatur för användande av slätväggig innerkanyl
- * Klar, mjuk, flexibel krage med storleksbeteckningar
- * Distalt fasad kanylände för slät övergång till obturator
- * Violet färgad låsbar obturator med innerlumen för ex ledare
- * Steril Etylenoxid - Pyrogenfri
- * Engångs
- * Latexfri
- * MR-kompatibel



Förpackningen innehåller;

- * Trakealkanyl
- * 2 st innerkanyler
- * Nackband
- * 1 st rengöringsborste till innerkanyl
- * Röntgentät

Benämning	Vårt Art nr	Enhet	Avdfp	Trpfp
Trakealkanyl, soft seal kuff i-kanyl 6,0 mm	100/810/060	st	1	10
Trakealkanyl, soft seal kuff i-kanyl 7,0 mm	100/810/070	st	1	10
Trakealkanyl, soft seal kuff i-kanyl 7,5 mm	100/810/075	st	1	10
Trakealkanyl, soft seal kuff i-kanyl 8,0 mm	100/810/080	st	1	10
Trakealkanyl, soft seal kuff i-kanyl 8,5 mm	100/810/085	st	1	10
Trakealkanyl, soft seal kuff i-kanyl 9,0 mm	100/810/090	st	1	10
Trakealkanyl, soft seal kuff i-kanyl 10 mm	100/810/100	st	1	10
Trakealkanyl m kuff o fen i-kanyl 6,0 mm	100/812/060	st	1	10
Trakealkanyl m kuff o fen i-kanyl 7,0 mm	100/812/070	st	1	10
Trakealkanyl m kuff o fen i-kanyl 7,5 mm	100/812/075	st	1	10
Trakealkanyl m kuff o fen i-kanyl 8,0 mm	100/812/080	st	1	10
Trakealkanyl m kuff o fen i-kanyl 8,5 mm	100/812/085	st	1	10
Trakealkanyl m kuff o fen i-kanyl 9,0 mm	100/812/090	st	1	10