

## PRODUKTINFORMATION- BIVONA®



### BIVONA® TRAKEALKANYL I SILIKON

- Det går att välja på "V" eller Rak krage
- Anpassad obturator med rundad spets
- En fast 15mm konnektor
- Mjuka, följsamma vingar med hål för fixering av kanylen med bomullsband eller nackband
- OD, ID & kanyllängd skrivet på kragen
- "SuperSlick" beläggning
- Ingen cuff, så inget skydd mot aspiration
- Knickresistent
- Armeringen innehåller inte järn, vilket ger mindre störningar vid MR undersökning

» CE märkt

» Material Silikon

» Steriliseringsmetod EO

» Enpatient

» MR-kompatibel

Benämning	ID (mm)	OD (mm)	Längd(mm)	Vårt Art nr	Enhet	Avdfp	Trpf
Bivona® trakealkanyl	2.5mm	4.0mm	38.0mm	60P025	st	1	1
Bivona® trakealkanyl	3.0mm	4.7mm	39.0mm	60P030	st	1	1
Bivona® trakealkanyl	3.5mm	5.3mm	40.0mm	60P035	st	1	1
Bivona® trakealkanyl	4.0mm	6.0mm	41.0mm	60P040	st	1	1
Bivona® trakealkanyl	4.5mm	6.7mm	42.0mm	60P045	st	1	1
Bivona® trakealkanyl	5.0mm	7.3mm	44.0mm	60P050	st	1	1
Bivona® trakealkanyl	5.5mm	8.0mm	46.0mm	60P055	st	1	1

PORTEX®

#### Smiths Medical Sverige AB

Box 1143, 164 22 Kista

Tel. 08- 594 772 50

Fax. 08-594 772 69

[www.smiths-medical.com](http://www.smiths-medical.com)

Vid frågor kontakta er  
lokala säljrepresentant

Smiths Medical is part of the global technology business Smiths Group plc. Product(s) described may not be licensed or available for sale in all countries. The Smiths Medical design mark is a trademark of Smiths Medical. The symbol © indicates the trademark is registered in the U.S. Patent and Trademark Office and certain other countries. ©2012 Smiths Medical. All rights reserved. 20121005SV-0001