

Gas sampling lines, standard

Vyaire gas sampling lines offer accurate CO₂ and anesthetic agent measurements and have been validated for use with the following GE Healthcare equipment: compact airway modules, CARESCAPE[™] respiratory modules and E-MiniC.



Cat. no.	Description	Qty.
5640I	Gas sampling line, 10', M/M	10
5640CEI	Co-extruded CO ₂ line male, 10'	20
5630I	Gas sampling line, 10', M/M with 0.8 filter	20
5640	Gas sampling line, 10', M/F	1
5630	0.8 micron filter, CO ₂ line, M/F connector	5
73325-HEL	Gas sampling line, 10', M/F	100
73319-HEL	Gas sampling line, 10', M/M	20
73318-HEL	Gas sampling line, 7', M/M	5
733199	Gas sampling line, 10', M/M	20
733188	Gas sampling line, 7', M/M	1
733170	Gas sampling line, 10', M/M with elbow	20
73306	Gas sampling line, 20', M/M	10
8004463	Exhaust line with colder	10
8004462	Exhaust line colder fit	100
M1031274	Gas sampling line, 7', M/F	10