

# VICRYL\*

(Polyglactin 910) Braided and Monofilament

## SUTURE

<b>Description:</b>	<b>Braided</b>	For easy handling and secure knot tying.
	<b>Coated</b>	For smooth passage through tissue and easy knot tie down.
	<b>Synthetic</b>	For minimal tissue reaction.
	<b>Absorbable</b>	By hydrolysis. Essentially complete in 56-70 days.

<b>Tensile strength:</b>	<b>Post implantation</b>	<b>Approximate % original strength remaining</b>
	14 days	75%
	21 days	50% (6-0 and larger) 40% (7-0 and smaller)
	28 days	25% (6-0 and larger)

**Color:** Violet or undyed.

**Range:** 10-0 to 2 (USP). Supplied as needed sutures and ligatures. Sizes 10-0 and 9-0 are monofilament in structure.












**Indications:** VICRYL\* sutures are intended for use in soft tissue approximation and/or ligation, including use in ophthalmic surgery, peripheral nerve adaptation and microsurgery for vessels less than 2 mm in diameter.

Typical areas of use include:







- Fascia closure - due to the unique 28-day claim
- Subcutaneous fat
- Joint capsule
- Uterus








Full details in the Instruction for Use included in every package.








1/2 Circle Taper Point Needles							
Needle	Gauge	Color	Length (cm)	Other information	Number per box	Code	
 TF-6 6,5 mm	8-0	●	30		12	V1027G	
	7-0	●	30		36	V1028H	
 TF-1 10 mm	5-0	●	70		36	V133H	
 RB-2 13 mm	6-0	●	70		36	PZ2571H	
 TF 13 mm	6-0	●	70		36	V926H	
	5-0	●	70		36	V925H	
	4-0	●	70		36	V924H	
 RB-1 17 mm	6-0	●	70		36	V302H	
	5-0	●	70		36	V303H	
	5-0	○	70		36	V213H	
	4-0	●	70		36	V304H	
	4-0	○	70		36	V214H	
	4-0	●	8x45	CONTROL RELEASE*	24	V734E	
	4-0	●	4x45	CONTROL RELEASE*	12	V735G	
	4-0	●	8x70	CONTROL RELEASE*	24	V799E	
	3-0	●	70		36	V305H	
	3-0	○	70		36	V215H	
	3-0	●	8x45	CONTROL RELEASE*	12	J713D	
2-0	●	70		36	V306H		
 RB-1 RB-1 17 mm	3-0	●	70		36	V667H	
 JRB-1  VISI-BLACK 17 mm	5-0	●	70	VISI-BLACK*	36	V3030H	
	4-0	●	70	VISI-BLACK*	36	V3040H	
 SH-1 22 mm	4-0	●	70		36	V310H	
	4-0	○	70		36	V426H	
	3-0	●	70		36	V311H	
	3-0	○	70		36	V219H	
	3-0	●	4x45	CONTROL RELEASE*	36	V782H	
	2-0	●	70		36	V312H	
	2-0	○	70		36	V220H	
	2-0	●	8x45	CONTROL RELEASE*	24	V792E	
 JB-1  VISI-BLACK 22 mm	4-0	●	70	VISI-BLACK*	36	V3100H	
	3-0	●	70	VISI-BLACK*	36	V3110H	
	2-0	●	70	VISI-BLACK*	36	V3120H	





**VICRYL\***  
(Polyglactin 910) Braided and Monofilament


 1/2 Circle Taper Point Needles						
Needle	Gauge	Color	Length (cm)	Other information	Number per box	Code
 SH 26 mm	4-0	●	70		36	V315H
	4-0	●	8x70	CONTROL RELEASE*	12	V783G
	3-0	●	70		36	V316H
	3-0	○	70		36	V416H
	3-0	●	90		36	V242H
	3-0	●	4x45	CONTROL RELEASE*	36	V780H
	3-0	●	8x45	CONTROL RELEASE*	12	J774D
	3-0	●	8x70	CONTROL RELEASE*	12	J784D
	2-0	●	70		36	V317H
	2-0	○	70		36	V224H
	2-0	●	90		36	V243H
	2-0	●	8x45	CONTROL RELEASE*	12	J775D
	0	●	70		36	V318H
	0	○	70		36	V225H
 JB VISI-BLACK 26 mm	4-0	●	70	VISI-BLACK*	36	V3150H
	3-0	●	70	VISI-BLACK*	36	V3160H
	3-0	○	70	VISI-BLACK*	36	V4416H
	2-0	●	70	VISI-BLACK*	36	V3170H
 MH-1 31 mm	3-0	●	70		36	V245H
	2-0	●	70		36	V246H
	0	●	70		36	V247H
 MH 36 mm	4-0	●	90		36	PG2592H
	3-0	●	70		36	V322H
	3-0	●	90		36	V100H
	2-0	●	70		36	V323H
	2-0	○	70		36	V226H
	2-0	●	90		36	V101H
	0	●	70		36	V324H
	0	●	90		36	V102H
	1	●	70		36	V325H
	1	○	70		36	V228H
 LH 40 mm	2-0	●	70		36	PH2574H





 1/2 Circle Taper Point Needles (Heavier Wire Diameter)						
Needle	Gauge	Color	Length (cm)	Other information	Number per box	Code
 UCL 13 mm	3-0	●	70		36	V113H
	0	●	70		36	V115H
 UCLX 16 mm	3-0	●	70		36	V13H
	2-0	●	70		36	V14H
	0	●	70		36	V15H
	1	●	70		36	V16H
 CT-3 22 mm	3-0	●	70		36	V327H
	2-0	●	70		36	V328H
	0	●	70		36	V329H
 CT-2 26 mm	3-0	●	70		36	V332H
	2-0	●	70		36	V333H
	2-0	○	70		36	V289H
	2-0	●	4x45	CONTROL RELEASE*	24	V7842E
	0	●	70		36	V334H
	0	○	70		36	V270H
	0	●	8x45	CONTROL RELEASE*	24	V744E
	0	●	4x70	CONTROL RELEASE*	24	4PK2575YE
	1	●	70		36	V335H
	2	●	70		36	V336H
 CT-1 36 mm	3-0	●	70		36	V338H
	3-0	○	90		36	V944H
	2-0	●	70		36	V339H
	2-0	●	90		36	V345H
	2-0	○	90		36	V945H
	2-0	○	8x70	CONTROL RELEASE*	12	J742D
	0	●	70		36	V340H
	0	●	90		36	V346H
	0	○	90		36	V946H
	0	●	8x45	CONTROL RELEASE*	12	J740D
	1	●	70		36	V341H
	1	●	90		36	V347H
	1	●	4x70	CONTROL RELEASE*	12	V736G
	2	●	90		36	V348H
2	○	90		36	V948H	
 MO-4 36 mm	0	●	70		36	V436H
	1	●	90		36	V255H





**VICRYL\***  
(Polyglactin 910) Braided and Monofilament



 <b>1/2 Circle Taper Point Needles (Heavier Wire Diameter)</b>							
Needle	Gauge	Color	Length (cm)	Other information	Number per box	Code	
 CT	40 mm	3-0	●	70		36	V350H
		2-0	●	70		36	V351H
		2-0	●	8x45	CONTROL RELEASE*	12	8PM2544YG
		0	●	70		36	V352H
		0	●	90		36	V358H
		0	●	150	LOOP	24	V1407E
		1	●	70		36	V353H
		1	●	90		36	V359H
		1	●	8x45	CONTROL RELEASE*	12	J753D
		2	●	90		36	V360H
		2	○	90		36	V960H
		 CTX	48 mm	3-0	●	70	
2-0	●			70		36	V363H
0	●			70		36	V364H
0	●			90		36	V370H
0	○			90		36	V170H
0	●			8x45	CONTROL RELEASE*	12	J764D
1	●			70		36	V365H
1	●			90		36	V371H
1	○			90		36	V366H
1	●			8x45	CONTROL RELEASE*	12	J765D
2	●			90		36	V372H
2	○			90		36	V367H
 TP-1	65 mm			2	●	4x70	
 XLH	70 mm	2-0	●	70		12	V581G
		0	●	70		12	V582G
		1	●	70		12	V583G


1/2 Circle ETHIGUARD* Blunt Point Needles						
Needle	Gauge	Color	Length (cm)	Other information	Number per box	Code
 MHB-1 31 mm	3-0	●	75		12	W9983
	2-0	●	75		12	W9984
 CTB-1 36 mm	2-0	●	75		12	W9985
	2-0	●	75	Heavy Needle	12	W9993
	0	●	75		12	W9986
	0	●	90	Heavy Needle	12	W9994
	1	●	90	Heavy Needle	12	W9995
 CTB 40 mm	0	●	75		12	W9988
	0	●	90	Heavy Needle	12	W9996
	1	●	90	Heavy Needle	12	W9997
 CTXB 48 mm	0	●	90		12	W9990
	1	●	90		12	W9991


3/8 Circle Blunt Point Needles						
Needle	Gauge	Color	Length (cm)	Other information	Number per box	Code
 BP-5 65 mm	1	●	100		12	W9391

5/8 Circle Taper Point Needles						
Needle	Gauge	Color	Length (cm)	Other information	Number per box	Code
 UR-6 26 mm	3-0	●	70		36	V601H
	2-0	●	70		36	V602H
	0	●	70		36	V603H
	1	●	70		36	VC2575H
 UR-6 26 mm	2-0	●	70		36	VCC2573H
 UR-5 36 mm	2-0	●	70		36	V375H
	0	●	70		36	V376H
	1	●	70		36	V377H
 UR-4 40 mm	2-0	●	70		36	V380H
	0	●	70		36	V381H

▼ 1/2 Circle TAPERCUT* Needles							
Needle	Gauge	Color	Length (cm)	Other information	Number per box	Code	
 V-37	40 mm	2-0	●	90		36	VS4020H
		0	●	90		36	VS4000H
 V-37	40 mm	1	●	70		36	ODD2574H
 V-40	48 mm	2-0	●	90		36	VS5020H
		0	●	70		36	VS5000H
		1	●	150	LOOP	24	V829E
 BPT-1	48 mm	6	●	75		12	V140G
		6	●	150		12	OH2528G



▼ 3/8 Circle TAPERCUT* Needles							
Needle	Gauge	Color	Length (cm)	Other information	Number per box	Code	
 V-18	13 mm	6-0	●	70		36	V991H
		5-0	●	70		36	V990H
 V-4	17 mm	4-0	○	70		12	V240G
		4-0	○	70		36	V240H




▼ 5/8 Circle TAPERCUT* Needles							
Needle	Gauge	Color	Length (cm)	Other information	Number per box	Code	
 UW-10	10 mm	6-0	●	30		12	W9172




▼ Compound Curve TAPERCUT* Needles							
Needle	Gauge	Color	Length (cm)	Other information	Number per box	Code	
 DA VISI-BLACK	25 mm	4-0	●	45	VISI-BLACK*	12	18531G

# VICRYL\*

(Polyglactin 910) Braided and Monofilament







● Fish Hook Needles						
Needle	Gauge	Color	Length (cm)	Other information	Number per box	Code
 PB-30 30 mm	2-0	●	70	Taper point	36	DF2570H
	1	●	70	Taper point	36	DF2573H
 V-38 31 mm	1	●	75	TAPERCUT*	12	W9335


● Straight Taper Point Needles						
Needle	Gauge	Color	Length (cm)	Other information	Number per box	Code
 ST-4 19 mm	4-0	●	70		36	BA2573H
	3-0	●	70		36	BA2574H
 ST-4 ST-4 19 mm	4-0	●	35		36	V620H
 ST 48 mm	3-0	●	70		36	BC2573H









▼ 1/2 Circle TAPERCUT* Needles						
Needle	Gauge	Color	Length (cm)	Other information	Number per box	Code
 V-5 17 mm	4-0	●	70		12	V994G
	4-0	●	70		36	V994H
	4-0	○	70		36	V230H
	3-0	●	70		36	V276H
 V-7 26 mm	3-0	●	70		12	V998G
	3-0	●	70		36	V998H
	2-0	●	70		36	V999H
	2-0	●	6x45	CONTROL RELEASE*	24	6OB2545YE
	0	●	70		36	V279H
 V-34 36 mm	3-0	●	90		36	V516H
	2-0	●	90		36	V517H
	0	●	90		36	V518H
	1	●	90		36	V519H
	1	○	90		36	V919H





**VICRYL\***  
(Polyglactin 910) Braided and Monofilament





▼ 1/2 Circle Reverse Cutting Needles (Heavier Wire Diameter)							
Needle	Gauge	Color	Length (cm)	Other information	Number per box	Code	
 X-1	22 mm	3-0	●	70		36	V460H
		3-0	○	70		36	ME2673H
		2-0	●	70		36	V461H
		2-0	○	75		36	JV459
 OS-4	22 mm	0	○	70	HEAVY NEEDLE	36	V694H
		2	●	75	HEAVY NEEDLE	12	W9287
 CP-2	26 mm	3-0	●	70		36	MF2572H
		2-0	○	70		36	V869H
		0	●	70		36	MF2574H
		0	●	4x45	CONTROL RELEASE	24	V870E
		1	○	70		36	MF2675H
 CP-1	36 mm	3-0	●	70		36	V465H
		2-0	●	70		36	V466H
		2-0	●	90		36	V472H
		2-0	○	90		36	JV1472
		0	●	70		36	V467H
		0	●	90		36	V473H
		0	●	4x45	CONTROL RELEASE	12	V475G
		1	●	70		36	V468H
		1	●	90		36	V474H
		2	●	90		36	V482H
 CP	40 mm	2-0	●	70		36	V478H
		0	●	70		36	V479H
		0	●	90		36	V485H
		1	●	70		36	V480H
		1	●	90		36	V486H
 CPX	48 mm	0	●	90		36	V1057H
		1	●	90		36	V1058H
		2	●	90		36	V1059H




▼ 1/2 Circle Reverse Cutting PRIME* Needles							
Needle	Gauge	Color	Length (cm)	Other information	Number per box	Code	
 PS-4	16 mm	5-0	○	45	Compound curve	12	W9436
		4-0	○	45	MULTIPASS*	36	MPV507H
		4-0	○	45		12	V507G
		4-0	○	45	Compound curve	12	W9437





▼ 3/8 Circle Reverse Cutting PRIME* Needles							
Needle	Gauge	Color	Length (cm)	Other information	Number per box	Code	
 P-1	11 mm	6-0	○	45	MULTIPASS*	36	MPV489H
		6-0	○	45		36	V489H
		5-0	○	45	MULTIPASS*	36	MPV490H
		5-0	○	45		36	V490H
 P-3	13 mm	6-0	○	45	MULTIPASS*	36	MPV492H
		6-0	○	45		36	V492H
		5-0	○	45	MULTIPASS*	36	MPV493H
		5-0	○	45		36	V493H
		4-0	○	45	MULTIPASS*	36	MPV494H
		4-0	○	45		36	V494H
 PS-3	16 mm	5-0	○	45	MULTIPASS*	36	MPV500H
		4-0	○	45		24	W9508T
 PS-2	19 mm	5-0	○	45	MULTIPASS*	36	MPV495H
		4-0	○	45	MULTIPASS*	36	MPV496H
		4-0	○	45		24	W9515T
		3-0	○	45	MULTIPASS*	36	MPV497H
		3-0	○	45		36	V497H
 PS-2 S	19 mm	4-0	○	45	Slimmed	36	V496H
 PS-1	24 mm	3-0	○	45	MULTIPASS*	36	MPV683H
		3-0	○	45		36	V683H
		2-0	○	45	MULTIPASS*	36	MPV684H
		2-0	○	45		36	V684H
 PS	26 mm	3-0	○	70	MULTIPASS*	36	MPV9526H
 PSL	30 mm	2-0	●	70		36	V686H



▼ 3/8 Circle Reverse Cutting Needles							
Needle	Gauge	Color	Length (cm)	Other information	Number per box	Code	
 C-3	13 mm	6-0	●	45		36	V384H
		5-0	●	45		36	V385H
		4-0	●	45		36	V386H
 FS-3	16 mm	6-0	●	45		36	V387H
		5-0	●	45		36	V388H



**VICRYL\***  
(Polyglactin 910) Braided and Monofilament


▼ 3/8 Circle Reverse Cutting Needles						
Needle	Gauge	Color	Length (cm)	Other information	Number per box	Code
 FS-2 19 mm	5-0	○	45		36	V291H
	5-0	●	45		36	V391H
	5-0	○	70		36	V421H
	4-0	○	45		36	V292ZH
	4-0	●	45		12	V392ZG
	3-0	○	45		36	V293H
	3-0	●	45		12	V393G
	3-0	●	45		36	V393H
	3-0	○	70		36	V423H
	3-0	●	70		36	V398H
 FS-2S 19 mm	4-0	○	45	Slimmed	36	V292H
	4-0	●	45	Slimmed	36	V392H
	4-0	○	70	Slimmed	36	V422H
	4-0	●	70	Slimmed	36	V397H
 FS-1 24 mm	4-0	○	70		36	V441H
	4-0	●	70		36	V451H
	3-0	○	70		36	V442H
	3-0	●	70		36	V452H
	2-0	○	70		36	V443H
	2-0	●	70		36	V453H
 FSL 30 mm	3-0	○	70		36	V588H
	3-0	●	70		36	V585H
	2-0	○	75		36	JV589
	2-0	●	70		36	V586H
	0	○	70		36	V520H
	0	●	70		36	V587H

▲ 3/8 Circle Conventional Cutting PRIME <sup>®</sup> Needles						
Needle	Gauge	Color	Length (cm)	Other information	Number per box	Code
 PC-3 16 mm	5-0	○	45	MULTIPASS <sup>®</sup>	36	MPV9505H
	5-0	○	45		24	W9505T
	4-0	○	45	MULTIPASS <sup>®</sup>	36	MPV9506H
 PC-5 19 mm	3-0	○	70	MULTIPASS <sup>®</sup>	36	MPV9571H
 PC-25 26 mm	4-0	○	70	MULTIPASS <sup>®</sup>	36	MPV9580H
	4-0	○	75		24	W9580T
	3-0	○	70	MULTIPASS <sup>®</sup>	36	MPV9581H
	3-0	○	75		24	W9581T

▲ Straight Cutting Needles							
Needle	Gauge	Color	Length (cm)	Other information	Number per box	Code	
 SC-1 13 mm	4-0	●	45		36	V383H	
 KS-25 25 mm	3-0	●	70		12	V2530G	
 TS 48 mm	3-0	○	70		36	AB2673H	
 KS 60 mm	4-0	○	70		36	V662H	
	3-0	●	70		36	V523H	
	3-0	○	70		36	V663H	
	2-0	●	70		36	V524H	
	2-0	○	70		36	V664H	




▼ 1/2 Circle MICROPOINT* Spatula							
Needle	Gauge	Color	Length (cm)	Other information	Number per box	Code	
 GS-14 GS-14 5,5 mm	10-0	●	30		12	V450G	
	9-0	●	30		12	V449G	
	8-0	●	20		12	V574G	
 TG175-8 7 mm	7-0	●	45		12	W9565	


▼ 3/8 Circle MICROPOINT* Spatula							
Needle	Gauge	Color	Length (cm)	Other information	Number per box	Code	
 GS-9 6,5 mm	8-0	●	30		12	V120G	
	7-0	●	30		12	V119G	
 GS-9 GS-9 6,5 mm	10-0	●	30	Monofilament	12	V440G	
	9-0	●	30		12	V549G	
	9-0	●	30	Monofilament	12	V439G	
	8-0	●	20		12	V547G	
	8-0	●	30		12	V548G	
	7-0	●	45		12	V546G	


▼ 1/4 Circle MICROPOINT* Spatula							
Needle	Gauge	Color	Length (cm)	Other information	Number per box	Code	
 TG100-8 6,5 mm	6-0	●	45		12	JV7544	


# VICRYL\*


(Polyglactin 910) Braided and Monofilament




▼ 3/8 Circle MICROPOINT* Reverse Cutting Needles						
Needle	Gauge	Color	Length (cm)	Other information	Number per box	Code
 G-6	8 mm	6-0	●	45		12 V551G
 G-1	11 mm	6-0	●	45		12 V554G
 G-3	13 mm	5-0	●	45		12 V557G






▼ 3/8 Circle CS ULTIMA* Spatula						
Needle	Gauge	Color	Length (cm)	Other information	Number per box	Code
 UCS-35	6,2 mm	10-0	●	15	MONOFILAMENT	12 V429G

▼ 1/8 Circle CS ULTIMA* Spatula						
Needle	Gauge	Color	Length (cm)	Other information	Number per box	Code
 UCS-26	6,2 mm	10-0	●	10	MONOFILAMENT	12 V865G

▼ 3/8 Circle Spatula Needles						
Needle	Gauge	Color	Length (cm)	Other information	Number per box	Code
 SPAT-10	11 mm	5-0	○	45		12 W9761

▼ 1/2 Circle Spatula Needles						
Needle	Gauge	Color	Length (cm)	Other information	Number per box	Code
 S-22	8 mm	6-0	●	45		12 W9566
		6-0	○	45		12 W9756

▼ 1/4 Circle Spatula Needles						
Needle	Gauge	Color	Length (cm)	Other information	Number per box	Code
 S-4	8 mm	4-0	●	45		12 V572G
 S-29	8 mm	6-0	●	45		12 V555G
 S-14	8 mm	6-0	●	45		12 V570G
		6-0	●	45	VISI-BLACK*	12 V575G
		6-0	○	45		12 V670G
		6-0	●	20		12 W9562
		5-0	●	45	VISI-BLACK*	12 V577G
		5-0	●	45		12 V571G
		5-0	○	45		12 V671G

 3/8 Circle Taper Point Needles							
Needle	Gauge	Color	Length (cm)	Other information	Number per box	Code	
BV-6 	4,8 mm	10-0	●	10	Monofilament	12	V970G
BV-4 	4,9 mm	10-0	●	15	Monofilament	12	V971G
BV130-5 	6,5 mm	9-0	●	15	Monofilament	12	V239G
		8-0	●	20		12	V400G
BV-1 	9,3 mm	6-0	●	45		12	V301G

Ligatures							
Needle	Gauge	Color	Length (cm)	Other information	Number per box	Code	
LIGAPAK*	4-0	●	250	Reel	12	V1204G	
	3-0	●	250	Reel	12	V1205G	
	2-0	●	250	Reel	12	V1206G	
	0	●	250	Reel	12	V1207G	
	1	●	250	Reel	12	V1208G	
	4-0	●	125	Reel	12	V204G	
	3-0	●	125	Reel	12	V205G	
	2-0	●	125	Reel	12	V206G	
	0	●	125	Reel	12	V207G	
SUTUPAK*	3-0	●	140		36	V614H	
	2-0	●	140		36	V615H	
	0	●	140		36	V616H	
	1	●	140		36	V617H	
	2	●	140		36	V618H	
SUTUPAK*	4-0	●	2x70		24	V623E	
	4-0	○	2x70		24	V651E	
	3-0	●	2x70		24	V624E	
	3-0	○	2x70		24	V652E	
	2-0	●	2x70		24	V625E	
	2-0	○	2x70		24	V653E	
	0	●	2x70		24	V626E	
	0	○	2x70		24	V654E	
	1	●	2x70		36	V627H	
	2	●	2x70		36	V628H	

**VICRYL\***  
(Polyglactin 910) Braided and Monofilament

Ligatures						
Needle	Gauge	Color	Length (cm)	Other information	Number per box	Code
SUTUPAK*	5-0	●	3x45		36	V632H
	5-0	○	3x45		36	V642H
	4-0	●	3x45		36	V633H
	4-0	○	3x45		36	V643H
	3-0	●	3x45		36	V634H
	3-0	○	3x45		36	V644H
	2-0	●	3x45		36	V635H
	2-0	○	3x45		36	V645H
	0	●	3x45		24	V636E
	1	●	3x45		24	V637E
	1	○	3x45		24	V647E
	2	●	3x45		36	V638H
SUTUPAK*	4-0	●	5x70		24	V1214E
	3-0	●	5x70		24	V1215E
	2-0	●	5x70		24	V1216E
	0	●	5x70		24	V1217E
	0	○	5x70		24	V965E
	1	●	5x70		36	V1218H
	2	●	5x70		36	V1219H
SUTUPAK*	5-0	●	6x45		36	V1223H
	4-0	●	6x45		36	V1224H
	3-0	●	6x45		36	V1225H
	3-0	○	6x45		36	V978H
	2-0	●	6x45		36	V1226H
	0	●	6x45		24	V1227E
	1	●	6x45		24	V1228E
	2	●	6x45		36	V1229H
SUTUPAK*	4-0	●	12x45		36	V903H
	3-0	●	12x45		36	V904H
	3-0	○	12x45		36	V910H
	2-0	●	12x45		24	V905E
	2-0	○	12x45		24	V911E
	0	●	12x45		24	V908E
	0	○	12x45		24	V912E
	1	●	12x45		24	V907E