


# Tubings


For hygienic transportation of secretions, reliable tubings complete the high quality system. Depending on the application, Medela offers various materials, diameters and connection pieces with the tubing.

## Silicone tubings



Ref. 077.0901 (1 m)  
 Ref. 077.0902 (2 m)  
 Ref. 077.0905 (5 m)  
 Ref. 077.0900 (25 m)

**Silicone tubings**  
 – Ø 7 x 12 mm  
 – Lengths 1, 2, 5, or 25 m



Ref. 077.0912 (60 cm)  
 Ref. 077.0913 (100 cm)

**Silicone tubings with 1 coupling piece**  
 – Ø 7 x 12 mm  
 – Lengths 60 or 100 cm



Ref. 077.0922 (60 cm)  
 Ref. 077.0921 (30 cm)

**Silicone tubings with 2 coupling pieces**  
 – Ø 7 x 12 mm  
 – Lengths 30 or 60 cm



Ref. 077.0931  
 Ref. 077.0185 (60 cm)

**Silicone tubing with 1 coupling piece and 1 coupling piece 90° angle**  
 – Ø 7 x 12 mm  
 – Length 25cm



Ref. 077.0970

**Silicone tubing**  
 – Ø 5 x 10 mm  
 – Length 25 m

## Disposable tubing



Ref. 077.0951

**Disposable tubing with fingertip**  
 – PVC, sterile  
 – Ø 6 x 9 mm  
 – Length 180 cm



Ref. 077.0965  
 Ref. 077.0952  
 Ref. 077.0191

**Disposable tubing with funnel connectors, Ø 5.5 x 8.5 mm**  
 – with two funnel connections, 150 cm  
 – with one funnel connection (non-sterile), 180 cm  
 – with one funnel connection, 300 cm



Ref. 077.0184

**Disposable tubing with universal connectors**  
 – for use with Single Use Bird Cup or reusable Bird Cups  
 – PVC, sterile  
 – Ø 5.5 x 8.3 mm  
 – Length 150 cm