

RT-series High Performance Breathing Circuits

Product Code: See table below

- » Superior condensate control in varying environments using spiral wire technology
- » Consistent compliance and low resistance minimize work of breathing
- » Probe latching designed to ensure correct orientation of probe for effective feedback
- » Circuit supplied preassembled
- » Lightweight and flexible



Adult breathing circuits

Part No.	Description	Quantity
RT200	Dual limb adult breathing circuit kit	10/box
RT204	Dual limb adult breathing circuit kit with pressure line	10/box
RT100	Dual limb adult breathing circuit (excl. chamber)	10/box
RT104	Dual limb adult breathing circuit (excl. chamber) with pressure line	10/box
RT205	Dual limb adult breathing circuit kit with expiratory water-trap	10/box
RT206	Dual limb adult breathing circuit kit with expiratory water-trap and pressure line	10/box
RT105	Dual limb adult breathing circuit (excl. chamber) with expiratory water-trap	10/box
RT106	Dual limb adult breathing circuit (excl. chamber) with expiratory water-trap and pressure line	10/box
RT319	Single limb adult breathing circuit kit with low resistance to flow inspiratory limb, flexible extension, exhalation port, inspiratory filter and pressure line	10/box
RT201	Single limb adult breathing circuit kit with Wye-piece and CPAP extension	10/box
RT101	Single limb adult breathing circuit (excl. chamber) with Wye-piece and CPAP extension	10/box
RT408	Single limb adult breathing circuit kit with flexible extension and Venturi air entrainer	10/box