

PROCARE®

California Tracheotomy Collar

Pre-formed molded foam tracheotomy collar provides support for both extrication and rehabilitation. Two-piece polyethylene foam lamination. Collar maintains head in neutral position and supports occipital and mandibular areas. Available in four (4) heights. Collar can be incorporated into a plaster body jacket.

Indication: Ideal for cervical spine immobilization and extraction.



Height*				
2.25"	3.25"	4.25"	5.25"	Size
79-83122				Ped./XS
79-83123	79-83133	79-83143	79-83153	S
79-83125	79-83135	79-83145	79-83155	M
79-83127	79-83137	79-83147	79-83157	L

XS = 7" - 10"; S = 10" - 13"; M = 13" - 16"; L = 16" - 19".

**Height = measurement from sternal notch to chin.*