



# LiNA Cool Brain

*Surgical headband with cooling effect*

## LiNA Cool Brain

The LiNA Cool Brain is a pre-moistened surgical headband with an automatic cooling effect to keep the surgeon and other OR staff cool and dry during surgery. Furthermore, the Cool Brain is soft and comfortable to wear and equipped with one-size elastic band to fit all head sizes.

The Cool Brain is made of PVA Foam which has an extremely large surface. The large surface allows the moisture to evaporate, giving the desired cooling effect. During surgery, the foam will absorb any perspiration, to continue the cooling.

The Cool Brain comes sterile and packed in a special aluminium pouch retaining moisture inside the headband for more than 3 years from date of production. The moisture consists of 99.8% water and 0.2% preservative.

- *Specially designed aluminium pouch retains moisture in headband even after 3 years*
- *Pre-moistened with one-size elastic band*
- *Ready for use*
- *Sterile*



### Ordering information:

LiNA Cool Brain

CE  
0459

REF.	Product	Packing
CB-200	<b>LiNA Cool Brain</b> - Surgical Headband made of PVA foam with pore size 130 micron. Pre-moistened and ready for use. Single use.	Box of 25, sterile

*LiNA distributor in your country:*

