

# Finn Chambers®

THE  
GOLD  
STANDARD  
for CONTACT ALLERGY  
PATCH TESTING



 **SmartPractice**

CALL +1 602.225.0595 | VISIT [smartpractice.com/dermatology](http://smartpractice.com/dermatology)

# Finn Chambers® The Gold Standard for Contact Allergy Patch Testing


## Finn Chambers® AQUA

- **Water Resistant!** Showers and exercise no longer interfere with adherence or results
  - **Comfortable!** Elastic polyurethane tape allows natural movement without pulling
  - **Convenient!** Pre-lined chambers allow for easy loading of allergens
- 100x10 patch test strips/box ..... Part # AL16630GEN



## Finn Chambers® on Scanpor®

- **Reliable!** The gold standard with verifiable occlusion and dependable adhesion
  - **Economical!** Flexible sizing and allergen-sparing chamber depth
  - **Patient Friendly!** Hypoallergenic Scanpor tape causes minimal skin irritation
  - **Aluminum** chambers are already mounted on tape for optimized skin adhesion, comfort, patient acceptability and ease of use
- 8 mm..... 25x10 patch test strips/box ..... Part # AL11111GEN  
 8 mm..... 100x10 patch test strips/box ..... Part # AL11133GEN  
 8 mm..... 50x1 patch test strips/box ..... Part # AL12121GEN  
 12 mm..... 200x1 patch test strips/box ..... Part # AL25150GEN
- Also available in **polypropylene-coated** chambers when using mercury solutions
- 8 mm..... 100x10 patch test strips/box ..... Part # AL11233GEN  
 8 mm..... 150x1 patch test strips/box ..... Part # AL14242GEN  
 12 mm..... 200x1 patch test strips/box ..... Part # AL25250GEN  
 18 mm..... 100x1 patch test strips/box ..... Part # AL35230GEN



## Finn Chambers® OPEN

- **Convenient!** Open design allows for easy allergen loading
  - **Reliable!** Dependable adhesion you expect from Finn Chambers
- 50 (5x2) patch test strips/box ..... Part # AL16121GEN

## Finn Chambers® ATOPY

- **Reliable!** Research has shown 12mm chambers are associated with a higher sensitivity (true positive) when used for atopic patch testing with food
- **Convenient!** 5 chambers on a sheet makes for easy loading and application of allergens for patch testing .....Part # AL23120GEN



## Loose Chambers

- 8mm .....1,000/bag ..... Part # AL4113GEN  
 12mm .....200/bag ..... Part # AL4212GEN  
 18mm .....100/bag ..... Part # AL4321GEN



EC REP

Dr. Ebeling & Assoc. GmbH  
 Isestrasse 5 • 20144 Hamburg • Germany  
**EC IMPORTER**  
 SmartPractice Europe GmbH  
 Boevemannstrasse 8 • 48268 Greven • Germany  
 +49 40 524707800

**SmartPractice**  
 3400 E McDowell Rd • Phoenix, AZ 85008 • U.S.A.

CALL +1 602.225.0595 | VISIT [smartpractice.com/dermatology](http://smartpractice.com/dermatology)