

ETHILON*

Monofilament Polyamide 6

SUTURE

Description: **Monofilament** For smooth passage through tissue.
Synthetic For minimal tissue reaction.
Non-absorbable Provides prolonged tensile strength retention in tissue.

Color: Blue, black or clear.








Range: 11-0 to 2 (USP). Supplied as needed sutures armed with, for example, MULTIPASS* and PRIME* needles. ETHILON* is also available as ligatures.

Indications: ETHILON* (black and clear) sutures are intended for use in general soft tissue approximation and/or ligation, including use in cardiovascular, ophthalmic, microsurgical and neurosurgical procedures. Blue ETHILON* sutures are intended for use in skin closure. Typical areas of use include:




- Cuticular sutures
- Nerve adaptation
- Ophthalmology









Full details in the Instruction for Use included in every package.












3/8 Circle Reverse Cutting PRIME* Needles							
Needle		Gauge	Suture Color	Length (cm)	Other information	Number per box	Code
 P-6	8 mm	7-0	●	45		36	1647H
		7-0	●	45		36	1696H
 P-1	11 mm	7-0	●	45	MULTIPASS*	36	MPE1696H
		6-0	●	45		36	697H
		6-0	●	45	MULTIPASS*	36	MPE697H
		6-0	●	45		36	1698H
 P-3	13 mm	6-0	●	45		36	1698H
		6-0	●	45	MULTIPASS*	36	MPE1698H
		5-0	●	45		36	698H
		5-0	●	45		36	EH7800H
		5-0	●	45	MULTIPASS*	36	MPE698H
		5-0	●	45	MULTIPASS*	36	MPE7800H
		5-0	●	45		24	W1852T
		4-0	○	45		36	691H
		4-0	●	45		36	699H
		4-0	●	45		36	EH7801H
		4-0	○	45	MULTIPASS*	36	MPE691H
		4-0	●	45	MULTIPASS*	36	MPE699H
		4-0	●	45	MULTIPASS*	36	MPE7801H
		 PS-3	16 mm	6-0	●	45	MULTIPASS*
6-0	●			45		24	W1614T
5-0	●			45		36	EH7930H
5-0	●			45	MULTIPASS*	36	MPE7930H
4-0	●			45		36	EH7931H
4-0	●			45	MULTIPASS*	36	MPE7931H
4-0	●			45		24	W1861T
 PS-2	19 mm	5-0	●	45		36	EH7804H
		5-0	●	45	MULTIPASS*	36	MPE7804H
		4-0	●	45		36	EH7162H
		4-0	●	45		24	W1620T
		4-0	●	45	MULTIPASS*	36	MPE1667H
		4-0	●	45	MULTIPASS*	36	MPE1855H
		3-0	●	45		36	EH7803H
		3-0	●	45	MULTIPASS*	36	MPE1621H
 PS-2 S	19 mm	4-0	●	45	SLIMMED	36	1667H
 PS-1	24 mm	4-0	●	45		36	EH7957H
		4-0	●	75		36	EH7956H
		4-0	●	75	MULTIPASS*	36	MPE7956H
		3-0	●	45		24	W1857T
		3-0	●	75		36	EH7954H
		3-0	●	75	MULTIPASS*	36	MPE7954H

ETHILON*
Monofilament Polyamide 6





3/8 Circle Reverse Cutting PRIME* Needles						
Needle	Gauge	Suture Color	Length (cm)	Other information	Number per box	Code
 PS 26 mm	3-0	●	45	MULTIPASS*	36	MPE1625H
	3-0	●	45		24	W1625T
	3-0	●	75		24	W1685T
 PSL 30 mm	3-0	●	75		36	EH7113H
 PSLX 36 mm	3-0	●	75		24	W1627T


3/8 Circle Reverse Cutting Needles						
Needle	Gauge	Suture Color	Length (cm)	Other information	Number per box	Code
 G-2 13 mm	6-0	●	45		36	667H
	5-0	●	45		36	EH7945H
 M-2 13 mm	5-0	●	45		36	644H
 FS-3 16 mm	6-0	●	35		12	W506
	6-0	●	45	VISI-BLACK*	36	EH7177H
	5-0	●	45		36	1628H
	5-0	●	45		36	EH7823H
	4-0	●	45		36	EH7824H
 FS-2 19 mm	5-0	●	45		36	EH7790H
	4-0	●	45		36	EH7144H
	4-0	●	75		36	EH7143H
	3-0	●	45		36	EH7792H
 FS-2S 19 mm	4-0	○	45	SLIMMED	36	611H
	4-0	●	45	SLIMMED	36	EH7791H
 FS-1 24 mm	5-0	●	45		36	EH7640H
	4-0	●	45		36	EH7825H
	4-0	●	45		36	S2450
	4-0	●	75		36	FJ5272H
	3-0	●	45		36	EH7794H
	3-0	●	75		36	EH7795H
	2-0	●	45		36	EH7826H
	2-0	●	45		36	EH7826H
 FS 26 mm	4-0	●	45		36	FK5242H
	3-0	●	45		36	ES163H
	2-0	●	45		36	EH7796H
 FSL 30 mm	4-0	●	75		36	1670H
	3-0	●	45		36	EH7936H
	3-0	●	75		36	EH7936H
	2-0	●	75		36	EH7798H


3/8 Circle Reverse Cutting Needles						
Needle	Gauge	Suture Color	Length (cm)	Other information	Number per box	Code
 FSLX 36 mm	3-0	●	75		36	1673H
	2-0	●	75		36	1674H
	0	●	75		36	EH7809H
 LR-40 60 mm	2-0	●	100		12	W558
 LR 76 mm	2	●	3x50		12	3460TG



3/8 Circle Conventional Cutting PRIME [®] Needles						
Needle	Gauge	Suture Color	Length (cm)	Other information	Number per box	Code
 PC-1 13 mm	6-0	●	45		36	JM5140H
 PC-3 16 mm	6-0	●	45	MULTIPASS [®]	36	MPE1853H
	6-0	●	45		24	W1815T
	6-0	●	45		24	W1853T
	5-0	●	45	MULTIPASS [®]	36	MPE1816H
	5-0	●	45		24	W1816T
	5-0	●	45		24	W1854T
 PC-5 19 mm	4-0	●	45	MULTIPASS [®]	36	MPE1856H
	4-0	●	45		24	W1856T
 PC-25 26 mm	6-0	●	45		24	W1859T
	5-0	●	45		24	W1858T
	5-0	●	45		12	W532
 KP-3 90 mm	2-0	●	100		12	W795
	1	●	100		12	W797
	1	●	2x50	With retention tubing	12	W2797
 KP-3 90 mm	1	●	75		12	W689


ETHILON*
Monofilament Polyamide 6

▼ 1/2 Circle Reverse Cutting Needles (Heavier Wire Diameter)							
Needle	Gauge	Suture Color	Length (cm)	Other information	Number per box	Code	
 M-1	13 mm	5-0	●	45		36	648H
 J-1	17 mm	5-0	●	75		36	F2423
 X-1	22 mm	3-0	●	45		36	EH7943H
 CPX	48 mm	1	●	100		12	W738

















▲ Straight Cutting Needles							
Needle	Gauge	Suture Color	Length (cm)	Other information	Number per box	Code	
 KS	60 mm	3-0	●	75		36	627H
		2-0	●	75		36	628H

◎ 1/2 Circle Taper Point Needles							
Needle	Gauge	Suture Color	Length (cm)	Other information	Number per box	Code	
 MH-1	31 mm	0	●	150	Loop	12	W739




◎ 1/2 Circle Taper Point Needles (Heavier Wire Diameter)							
Needle	Gauge	Suture Color	Length (cm)	Other information	Number per box	Code	
 CT	40 mm	0	●	200	Loop	12	W760
 CTX	48 mm	1	●	150	Loop	12	W748

◎ Straight Taper Point Needles							
Needle	Gauge	Suture Color	Length (cm)	Other information	Number per box	Code	
 STAW	65 mm	2-0	●	75		24	EH7625E




 3/8 Circle Taper Point Needles


Needle	Gauge	Suture Color	Length (cm)	Other information	Number per box	Code	
 BV 140-3	3,7 mm	10-0	●	15	Non-sterile	12	NS2820
 BV130-3	3,7 mm	9-0	●	13		12	FG2819
 BV75-3	3,8 mm	10-0	●	13		12	W2870
		9-0	●	13		12	W2871
 BV50-3	3,8 mm	11-0	●	13		12	W2881
 BV130-4	4,7 mm	9-0	●	13		12	W2813
		8-0	●	13		12	W2812
 BVH100-3	5 mm	10-0	●	13		12	W2880
 BV75-4	5,1 mm	10-0	●	13		12	W2850
 BV-6	5,1 mm	11-0	●	10		12	EH7438G
 BV-4	5,1 mm	10-0	●	15		12	2814G
 BV100-4	5,1 mm	10-0	●	13		12	W2830
		9-0	●	13		12	W2829
		9-0	●	15		12	EH7448G
 BV-2	6,5 mm	10-0	●	15		12	2810G
		9-0	●	15		12	2809G
		8-0	●	15		12	2808G
 BV130-5	6,5 mm	10-0	●	13		12	W2810
		8-0	●	13		12	W2808
 BV130-5	6,5 mm	10-0	●	30		12	W2801
 BV175-6	8 mm	8-0	●	45		12	W8170
 BV-1	9,3 mm	7-0	●	45		36	EH7448H
 C-1	13 mm	6-0	●	75		36	EH7500H



 3/8 Circle CS ULTIMA* Spatula


Needle	Gauge	Suture Color	Length (cm)	Other information	Number per box	Code	
 UCS-35	6,2 mm	10-0	●	15		12	7060G
 CS-140-6	6,5 mm	10-0	●	30		12	U7003
 CS-140-6	6,5 mm	10-0	●	30		12	U7061


ETHILON*
Monofilament Polyamide 6


▼ 3/8 Circle MICROPOINT* Spatula						
Needle	Gauge	Suture Color	Length (cm)	Other information	Number per box	Code
CS-35-C 	6 mm	10-0	●	15	12	EH7980G
GS-9 	6,5 mm	10-0	●	30	12	1718G
		9-0	●	30	12	1717G
TG140-6  TG140-6	6,5 mm	11-0	●	30	12	W1780
		10-0	●	30	12	W1770
		9-0	●	30	12	W1769



▼ 1/2 Circle MICROPOINT* Reverse Cutting Needles						
Needle	Gauge	Suture Color	Length (cm)	Other information	Number per box	Code
G-7 	8 mm	8-0	●	45	12	W1765


▼ 1/2 Circle MICROPOINT* Spatula						
Needle	Gauge	Suture Color	Length (cm)	Other information	Number per box	Code
TG160-6 	5,5 mm	10-0	●	30	12	W1758
TG175-8 	7 mm	8-0	●	30	12	1716G

▼ 1/2 Circle CS ULTIMA* Spatula						
Needle	Gauge	Suture Color	Length (cm)	Other information	Number per box	Code
CS-160-6 	5,5 mm	10-0	●	30	12	U7000

▼ 1/4 Circle MICROPOINT* Spatula						
Needle	Gauge	Suture Color	Length (cm)	Other information	Number per box	Code
TG100-8 	6,5 mm	8-0	●	30	12	B7714

▼ 1/4 Circle Spatula Needles						
Needle	Gauge	Suture Color	Length (cm)	Other information	Number per box	Code
S-14 	8 mm	5-0	●	20	36	MN5132H

1/4 Circle TAPERCUT* Needles						
Needle	Gauge	Suture Color	Length (cm)	Other information	Number per box	Code
V100-3 	3,5 mm	10-0	●	13		W2899
		9-0	●	13		W2898
V150-4 	4,5 mm	10-0	●	15		W1701

Compound Curve CS ULTIMA* Spatula						
Needle	Gauge	Suture Color	Length (cm)	Other information	Number per box	Code
CS-C-6 	5 mm	10-0	●	30		U7009

Ligatures						
Needle	Gauge	Suture Color	Length (cm)	Other information	Number per box	Code
SUTUPAK*	4-0	●	3x45		36	EH6751H
		●	10x45		36	EH6562H
LIGAPAK*	2-0	●	250	Reel	20	F4354