

C-147 Hållare för 1000ml flaska (rund)

Informations blad för Handskhållare C-147

Artikel nr: C-147

Lätt och smidig hållare till 1 liters flaska som passar i vårt fäste även på våra linjalser/dörrstationer/iva och EU. Ändras förpackningen kan smidigt enbart hållaren bytas ut och fästet sitta kvar på väggen utan att behöva ta ner och skruvas bort och måla.

Man sätter snabbt och enkelt i flaskan!

Man byter enkelt hållare om man vill, Passa alla Sanaquas hållare/Fästen.

Komplett förpackning innehåller:

- väggfäste med dubbelhäftande tape
- rengöringsservett
- skruv
- plugg
- stålkorg





"Effektivt och hygieniskt - allt på sin plats vid varje tillfälle"

Our Multi Holders

Mounting rulers for Sanaqua wirebaskets and accessories. For horizontal mounting. Screw and plug included.



C-175-1

Ruler 2 holder
Length = 230 mm



C-175-2

Ruler long 2 holder
Length = 340 mm



C-174-1

Ruler 3 holder
Length = 342 mm



Mounting rulers for Sanaqua baskets and accessories. For vertical mounting. Screw and plug included.



C-165-1

Ruler 2 holder
Height = 301 mm



C-164-1

Ruler 3 holder
Height = 458 mm



C-162-1

Ruler 4 holder
Height = 615 mm



C-159-1

Ruler 5 holder
Height = 772 mm

