


<b>Product Name</b>	Scavenging connector, 22M - 30F	
<b>Reference numbers</b>	1971000	
<b>Product Type</b>	Connector	
<b>Manufacturer</b>	Intersurgical	
<b>Materials</b>	SB	
<b>Packaging</b>		
<b>Unit Packaging</b>	Individually packaged	
<b>Box quantity</b>	25	
<b>Pallet quantity</b>	6400	
<b>Instructions for Use</b>	Instructions for use with product	N/A
	Instructions for use with each box	N/A
<b>CE Class</b>	IIa	
<b>Applicable Standards</b>	BS EN ISO 13485: Medical devices — Quality management systems — Requirements for regulatory purposes	
	ISO 9001: Quality Management Systems - Requirements	
	ISO 15223-1: Medical Devices - Symbols to be Used with Medical Device Labels, Labelling, & Information to be Supplied	
	BS EN ISO 5356-1: Anaesthetic and respiratory equipment. Conical connectors. Cones and sockets	
<b>Sterile</b>	No	
<b>Shelf life</b>	5 Years	
<b>Disposal Methods</b>	As hospital guidelines	
<b>Re-use</b>		Single use only

## Product characteristics

### Technical information

Weight (g)	8.5
Internal volume (ml)	24.4



**Material  
Abbreviations**

ABS	Acrylonitrile Butadiene Styrene
AL	Aluminium
CR	Polychloroprene Rubber
HDPE	High Density Polyethylene
HIPS	High Impact Polystyrene
LDPE	Low Density Polyethylene
MBS	Methacrylate Butadiene Styrene
NRL	Natural Rubber Latex
PC	Polycarbonate
PE	Polyethylene
PP	Polypropylene
PS	Polystyrene
PSU	Polysulfone
PVC	Poly Vinyl Chloride
SB	Styrene Butadiene
Si	Silicone
SS	Stainless Steel
TPE	Thermoplastic Elastomer
SBC	Styrene Butadiene Copolymer

