

### Vital Signs laryngoscope handles

The Vital Signs reusable laryngoscope handle is ergonomic with a comfortable over-mold and lightweight design to fit the clinician's hand. It contains a durable, slip-resistant grip and a wide, sturdy base. Contains a bright LED light and no bulb replacement required. 5-year warranty on all handles. Efficient technology extends battery life beyond that of a standard incandescent bulb light source. Proprietary connection to Vital Signs blades. Batteries are not included.



Cat. no.	Description	Qty.
4556EU	Standard	1
4557VS	Stubby	1

### Vital Signs laryngoscope plastic blades

Our laryngoscope system is exclusively provided by Vyair with a proprietary secure connection. The plastic laryngoscope blades provide a cost-effective option to metal without compromising strength. These single-patient-use blades can be incinerated with no negative byproducts.



Cat. no.	Description	Qty.
4002VS	Macintosh 2	20
4003VS	Macintosh 3	20
4004	Macintosh 4	20
4510	Miller 0	20
4511	Miller 1	20
4512VS	Miller 2	20
4513VS	Miller 3	20
<b>High visibility</b>		
4502	Macintosh 2	20
4503	Macintosh 3	20
4504EU	Macintosh 4	20
<b>With sheath</b>		
4002S	Macintosh 2	20
4003S	Macintosh 3	20
4004S	Macintosh 4	20