

HANDHELD INTERMITTENT CATHETERS


LATEX
FREE



The Handheld Intermittent Catheter is fully adjustable for comfort and ease of use for self catheterisation.

Movable handle to ascertain the ideal fitting position.

The catheter has nelaton tip and made from PVC to minimise allergic reactions.



- Individually supplied sterile in easy to peel pouches
- Latex-Free
- CE Marked

IMP CODE	DESCRIPTION	CATHETER SIZE	HANDLE SIZE	NHS CODE	PACK QTY
INS/100/8	Intermittent Catheter	180mm x 8ch	115mm Length	-	100
INS/100/10	Intermittent Catheter	180mm x 10ch	115mm Length	-	100
INS/100/12	Intermittent Catheter	180mm x 12ch	115mm Length	-	100
INS/100/14	Intermittent Catheter	180mm x 14ch	115mm Length	-	100



INSIGHT
Medical Products Ltd

Insight Medical Products Ltd
Unit 3, Priory Industrial Estate,
Tetbury, Gloucestershire. GL8 8HZ

Tel: +44 (0)1666 500055
info@insightmedical.net
www.insightmedical.net

