


<b>Product Name</b>	Flexible straight catheter mount 22F-22M/15F, 170mm	
<b>Reference numbers</b>	3511000	
<b>Product Type</b>	Catheter Mount	
<b>Manufacturer</b>	Intersurgical	
<b>Materials</b>	SB	LDPE
<b>Packaging</b>		
<b>Unit Packaging</b>	Individually packaged	
<b>Box quantity</b>	50	
<b>Pallet quantity</b>	3200	
<b>Instructions for Use</b>	Instructions for use with product NO Instructions for use with each box NO	
<b>Certification</b>	CE0120	Class 2A
<b>Applicable Standards</b>	BS EN ISO 13485: Medical devices — Quality management systems — Requirements for regulatory purposes ISO 9001: Quality Management Systems - Requirements ISO 15223-1: Medical Devices - Symbols to be Used with Medical Device Labels, Labelling, & Information to be Supplied BS EN ISO 5356-1: Anaesthetic and respiratory equipment Conical connectors. Cones and sockets. BS 7143: Catheter mounts (flexible adaptors) for use with medical breathing systems	
<b>Sterile</b>	No	
<b>Shelf life</b>	5 Years	
<b>Disposal Methods</b>	As hospital guidelines	
<b>Re-use</b>	 Single use only	

## Product characteristics

### Technical information

Length (mm)	170
Weight (g)	15.8
Resistance to flow at 30 L/min (mbar)	0.2
Leakage at 60 mbar (ml/min)	0.0
Compressible volume (ml)	41.9
Connectors (Machine / Patient)	22F / 22M/15F



**Material  
Abbreviations**

ABS	Acrylonitrile Butadiene Styrene
AL	Aluminium
HDPE	High Density Polyethylene
HIPS	High Impact Polystyrene
LDPE	Low Density Polyethylene
MBS	Methacrylate Butadiene Styrene
NRL	Natural Rubber Latex
PC	Polycarbonate
PE	Polyethylene
PP	Polypropylene
PVC	Poly Vinyl Chloride
SB	Styrene Butadiene
Si	Silicone
SS	Stainless Steel
TPE	Thermoplastic Elastomer
SBC	Styrene Butadiene Copolymer

