



Appendix 1: Infusion Sets for Alaris™ Pumps (GP, VP, CC, GW/GW800 and SE, IVAC 590 series)

All products listed in Appendix 1, manufactured by BD Switzerland Sàrl, Terre Bonne Park – A4, Route de Crassier 17, Eysins, Switzerland unless otherwise indicated.

Product Catalogue (REFs)	Product Name
Infusion Sets for Alaris™ CC Syringe Pump	
G30402M	Alaris™ Products. Alaris™ CC Syringe Pump, Asena™ CC Syringe Pump, IVAC™ 701/711/770/P7000 Pumps. EXTENSION SET.
G30453V	Alaris™ Products. Alaris™ CC Syringe Pump, Asena™ CC Syringe Pump, IVAC™ P7000 Pumps. EXTENSION SET. LIGHT RESISTANT. LOW SORBING (PE/PVC)
G30653V	Alaris™ Products. Alaris™ CC Syringe Pump, Asena™ CC Syringe Pump, IVAC™ P7000 Pumps. EXTENSION SET. LIGHT RESISTANT.
Infusion Sets for Alaris™ GP Volumetric Pump	
60593	Alaris™ Products. Alaris™ GP Volumetric Pump. Alaris™ Safety Clamp. INFUSION SET. 15 µm FILTER.
60643	Alaris™ Products. Alaris™ GP Volumetric Pump. Alaris™ Safety Clamp. INFUSION SET. LIGHT RESISTANT. 15 µm FILTER.
60693	Alaris™ Products. Alaris™ GP Volumetric Pump. Alaris™ Safety Clamp. INFUSION SET. 15 µm FILTER.
60793	Alaris™ Products. Alaris™ GP Volumetric Pump. Alaris™ Safety Clamp. INFUSION SET. 15 µm FILTER.
60894	Alaris™ Products. Alaris™ GP Volumetric Pump. Alaris™ Safety Clamp. TRANSFUSION/BLOOD SET. 200 µm FILTER.
60895	Alaris™ Products. Alaris™ GP Volumetric Pump. Alaris™ Safety Clamp. TRANSFUSION/BLOOD SET. 200 µm FILTER.
60903	Alaris™ Products. Alaris™ GP Volumetric Pump. Alaris™ Safety Clamp. INFUSION SET. 15 µm FILTER.
60980	Alaris™ Products. Alaris™ GP Volumetric Pump. Alaris™ Safety Clamp. TRANSFUSION/BLOOD SET. 2 SPIKE. 200 µm FILTER.
60093E	Alaris™ Products. Alaris™ GP Volumetric Pump. Alaris™ Safety Clamp. SmartSite™ Needle-free System. INFUSION SET. 15 µm FILTER.
60123E	Alaris™ Products. Alaris™ GP Volumetric Pump. Alaris™ Safety Clamp. SmartSite™ Needle-free System. INFUSION SET. 1.2 µm FILTER.
60173E	Alaris™ Products. Alaris™ GP Volumetric Pump. Alaris™ Safety Clamp. SmartSite™ Needle-free System. INFUSION SET. NO FILTER.
60393E	Alaris™ Products. Alaris™ GP Volumetric Pump. Alaris™ Safety Clamp. SmartSite™ Needle-free System. TRANSFUSION / BLOOD SET. 200 µm FILTER.
60693E	Alaris™ Products. Alaris™ GP Volumetric Pump. Alaris™ Safety Clamp. SmartSite™ Needle-free System. INFUSION SET. 15 µm FILTER.
60793E	Alaris™ Products. Alaris™ GP Volumetric Pump. Alaris™ Safety Clamp. SmartSite™ Needle-free System. INFUSION SET. 15 µm FILTER.
60950E	Alaris™ Products. Alaris™ GP Volumetric Pump. Alaris™ Safety Clamp. SmartSite™ Needle-free System. MULTI-INFUSION SET. ONCOLOGY. VENTED. 15 µm FILTER.
60951E	Alaris™ Products. Alaris™ GP Volumetric Pump. Alaris™ Safety Clamp. SmartSite™ Needle-free System. MULTI-INFUSION SET. ONCOLOGY. VENTED. 15 µm FILTER.
60952E	Alaris™ Products. Alaris™ GP Volumetric Pump. Alaris™ Safety Clamp. SmartSite™ Needle-free System. MULTI-INFUSION SET. ONCOLOGY. LIGHT RESISTANT. VENTED. 15 µm FILTER.
60953V	Alaris™ Products. Alaris™ GP Volumetric Pump. Alaris™ Safety Clamp. INFUSION SET. LOW SORBING. 15 µm FILTER.
63200NYB	Alaris™ Products. Alaris™ GP Volumetric Pump. Alaris™ Safety Clamp. INFUSION SET. NO FILTER.
63401EB	Alaris™ Products. Alaris™ GP Volumetric Pump. Alaris™ Safety Clamp. SmartSite™ Needle-free System. INFUSION SET. NO FILTER.
63402BE	Alaris™ Products. Alaris™ GP Volumetric Pump. Alaris™ Safety Clamp. SmartSite™ Needle-free System. INFUSION SET. NO FILTER.
63420EB	Alaris™ Products. Alaris™ GP Volumetric Pump. Alaris™ Safety Clamp. SmartSite™ Needle-free System. INFUSION SET. NO FILTER.
63423BE	Alaris™ Products. Alaris™ GP Volumetric Pump. Alaris™ Safety Clamp. SmartSite™ Needle-free System. INFUSION SET. NO FILTER.
63477EB	Alaris™ Products. Alaris™ GP Volumetric Pump. Alaris™ Safety Clamp. SmartSite™ Needle-free System. TRANSFUSION/BLOOD SET. 2 SPIKE. 200 µm FILTER.
Infusion Sets for Alaris™ GW Volumetric Pump	