

Description

Screw cap container

Art. No.	75.562.300	No image available Keine Abbildung vorhanden
Volume	100 ml	
Diameter	57 mm	
Length	76 mm	
Base shape	flat base	
Type of closure	screw cap	
Material	polypropylene (PP)	
Colour	transparent	
Max. centrifugation speed	No details available	
Please refer to our notes on the product characteristics.		
pcs. / sales unit	300	
pcs. / sub-package	50	
<ul style="list-style-type: none"> ■ Non-sterile tube ■ Autoclavable to 121°C. ■ Screw cap enclosed. ■ Moulded graduation. 		