

# Micro Tubes with Attached Loop Screw Cap

For transport, storage and sample processing in the laboratory



- Attached loop screw cap avoids mix-ups, screw cap with O-ring
- Tube and cap from high-quality polypropylene, autoclavable at 121°C
- Knurls enable single-handed operation in special Sarstedt racks for micro tubes (racks on pp. 113 -115)

## Micro tubes 0.5 ml with skirted base, neutral cap

Order Number	Description	Packaging Pcs./Case
72.730.100	attached screw cap, non-assembled	5,000
72.730.107	with writing space, attached screw cap, non-assembled	5,000
72.730.105	screw cap assembled, sterile	1,000
72.730.106	with writing space, screw cap assembled, sterile	1,000



## Micro tubes 1.5 ml, neutral cap


Order Number	Description	Packaging Pcs./Case
72.692.100	attached screw cap, non-assembled	5,000
72.692.105	screw cap assembled, sterile	1,000

Storage boxes and special racks for micro tubes on p. 105 pp. 113-115.



## Micro tubes 2.0 ml, neutral cap

Order Number	Description	Packaging Pcs./Case
72.693.100	attached screw cap, non-assembled	5,000
72.693.105	screw cap assembled, sterile	1,000



## Micro tubes 2.0 ml with skirted base, cap neutral

Order Number	Description	Packaging Pcs./Case
72.694.100	attached screw cap, non-assembled	5,000
72.694.107	with writing space and graduations*, attached screw cap, non-assembled	5,000
72.694.105	screw cap assembled, sterile	1,000
72.694.106	with writing space and graduations*, attached screw cap, non-assembled	1,000

## Screw cap with separate retaining loop, neutral

Order Number: 65.712, packaging 500 Pcs./Bag, 5,000 Pcs./Case

## Colour coded inserts

(only fit attached caps)

Packaging 500 Pcs./Bag, 5,000 Pcs./Case

Order Number	Colour	Order Number	Colour
65.713	neutral	65.713.006	blue
65.713.002	red	65.713.007	violet
65.713.004	yellow	65.713.009	white
65.713.005	green	65.713.010	black



Length incl. cap – 49 mm, Ø below collar – 10.8 mm

\*For technical reasons, the graduations can only be printed up to 1.25 ml.