




Medline International France SAS

5 rue Charles Lindbergh
44110 Châteaubriant

www.medline.com/uk



Technical sheet

Surgical gown with fabric reinforcement

Reference(s)

9451CEA, 9508CEA, 9509CEA, 9512CEA, 9539CEA, 9571CEA, 9581CEA

Description

A sterile wrap 60x60 cm contains:

- 1 Optima® surgical gown, with fabric reinforcement, Velcro® adjustable neck
- Hand towels, 39x30 cm unfolded

Reference	Length (in cm)	Size	Qty of towels	Qty of gown(s) per pack factor	Unit/case
9451CEA	150	XLL	4	1	20
9508CEA	120	L	2	1	24
9509CEA	130	XL	2	1	22
9512CEA	150	XLL	2	1	22
9539CEA	130	XL	6	3	7
9571CEA	130	XXL	2	1	20
9581CEA	170	2XLL	2	1	19

Main materials

Optima® (polyester + woodpulp)

Reinforcement:

- Sleeves: Optima® (polyester + woodpulp)
- Chest: Wetlaid (cellulose)

Packaging

Peel pouch

2 polyliners contain all the units in a carton box

Manufacturer

Medline International France SAS

EN ISO 13485

Regulation

In compliance with Medical Devices Directive 93/42/EEC

MDD Class: I sterile



Latex free



0123



STERILE | EO

© 2011 Medline Industries Inc. or one of its subsidiaries.
All rights reserved.

LIT-95010CE-EN02

01.05.SUB.04