

ABG SYRINGES

Product Number: 364314



MSDS

Certificate of Analysis

364314 - 3 mL BD Preset™ syringe with BD Luer-Lok™ tip. 22 G x 1 in. needle. Green plunger stopper, tip cap, and blue cube. CE marked. (100/ca)

Key Product Features	
Availability	Canada
Product type	Blood gas analysis syringe
Brand	BD Preset™, BD Luer-Lok™
Closure color	Green
DNR (Latex) free?	Yes
Draw volume	3 mL
Needle gauge	22 G
Needle length	1.00 in.
Sampling technique	BD Preset™
Tip	Luer-Lok
Type	BD Preset™