




Medline International France SAS

5 rue Charles Lindbergh
44110 Châteaubriant

www.medline.com/uk



Technical sheet

Extremity drape

Reference(s)

29414CE

Description

A sterile wrap 76x76 cm contains:

- 1 Tiburon® extremity drape, 221x326 cm, 74x114 cm reinforcement, 6 cm diameter elastic fenestration, 6 cord holding tabs

Main materials

Tiburon® (polypropylene + polyethylene)

Polyrayon reinforcement 64g/m² (rayon + polyester)

Pack factor

12 units per carton

Packaging

Vented sterilization pouch

1 polyliner contains all the units in a carton box

Manufacturer

Medline International France SAS

EN ISO 13485

Regulation

In compliance with Medical Devices Directive 93/42/EEC

MDD Class: I

HP: High performance according to EN 13795 standard



Latex free



0123



STERILE EO



© 2011 Medline Industries Inc. or one of its subsidiaries.
All rights reserved.

LIT-29414CE-EN01

01.05.SUB.04