

Alaris™ GP Series Primary infusion set with 15µm chamber filter



Alaris™ GP Series Primary infusion set with 15µm chamber filter

60593

- Set Length 265 cm
- Priming Volume 24 ml
- Bore 3 mm
- 1 Spiked Drip Chamber(s)
(Vented/Unvented)
- Drops per ml 20
- Hydrophobic Air Inlet Filter (with Cap)
- 15 µm Drip Chamber Filter
- 1 Roller Clamp(s)
- Alaris™ Safety Clamp
- Silicone pump segment
- 1 Male Luer(s)
- Polyvinyl Chloride (PVC)
- Natural rubber latex is not part of the material formulation
- DEHP is not part of the material formulation
- Quantity 80



1 Male Luer(s)



86