

Alaris™ GP Plus Series Blood Set with 200µm chamber filter



Alaris™ GP Plus Series Blood Set with 200µm chamber filter

63477EB

- Needle-Free Connectors
- Set Length 300 cm
- Priming Volume 32 ml
- Bore 3 mm
- 2 Spike(s) (Unvented)
- In-Line Drip Chamber
- Drops per ml 20
- 200 µm Drip Chamber Filter
- 1 SmartSite® Valve Y Port
- 1 Roller Clamp(s)
- 2 Slide Clamp(s)
- Alaris™ Safety Clamp
- Pump Segment
- 1 Male Luer(s)
- Natural rubber latex is not part of the material formulation
- DEHP is not part of the material formulation
- Quantity 20



1 SmartSite® Valve Y Port 1 Male Luer(s)

0086