

Portex Tracheostomy Tubes

Portex® Blue Line® Cuffed Adjustable Flange Tracheostomy Tubes

Not available in USA

[INQUIRE ABOUT THIS PRODUCT](#)

PORTEX

The tube features a unique flange which can be adjusted by a simple locking device. Designed for use in patients with abnormally large anatomies, such as glottic oedema or obesity, who require an extended tracheostomy tube from flange to bend.

The device can be adjusted along the tube to obtain the most comfortable and safe tube position.

Product Features

- ▶ 15mm Connector complies with ISO 5356-1. Allows connection to breathing system or a Heat & Moisture Exchange, Thermovent T (100/570)
- ▶ Smooth obturator facilitates insertion of tracheostomy tube
- ▶ Thermosensitive PVC Flange
- ▶ Anatomically shaped tube
- ▶ Includes cotton neck tapes



Ordering Information:

Box quantity: 2



[Symbol Key](#)

[Product Codes](#)

[Documents](#)

[Associated Products](#)

Product Specifications



Cuffed Blue Line® Adjustable Flange Tracheostomy Tubes

Reorder Code	ID (mm)	OD (mm)	Uncuffed Code
100/523/060	6.0mm	8.2mm	100/526/060
100/523/070	7.0mm	9.6mm	100/526/070
100/523/080	8.0mm	10.9mm	100/526/080
100/523/090	9.0mm	12.3mm	100/526/090
100/523/100	10.0mm	13.7mm	100/526/100

Notes: No Inner cannulae available for these tubes