

NEW! 8.3F KITS



The Right Combination

*Everything necessary for the procedure is on
your shelf in an organized, ready-to-use tray.*

PERICARDIOCENTESIS KIT™

BE PREPARED!

Merit's Pericardiocentesis Kits offer:

THE RIGHT COMBINATION

Everything necessary for the procedure is on your shelf in an organized, ready-to-use tray.

TIME SAVINGS

Time previously spent assembling a kit from scratch can now be focused toward the patient.

COST SAVINGS

Ordering the ready-to-use kit saves on personnel time and lowers overall component cost.

PROCEDURE SPECIFIC DESIGN

High-flow catheters are developed with pericardiocentesis procedures in mind.

ORDERING INFORMATION

Cat. No.	Description
PC801	8.3F Pigtail Catheter
PC802	8.3F Straight Catheter
PC803EU	8.3F Pigtail Catheter, Europe (CE Mark Pending)
PC804EU	8.3F Straight Catheter, Europe (CE Mark Pending)
PC100	6F Pigtail Catheter, Europe (CE Mark)
PC200	6F Straight Catheter, Europe (CE Mark)
PC101	6F Pigtail Catheter
PC201	6F Straight Catheter



U.S. Patent No. 5,403,337 Other U.S. and Foreign Patents Pending



Merit Medical Systems, Inc. • 1600 West Merit Parkway, • South Jordan, Utah 84095 • 1-801-253-1600 • 1-800-35-MERIT
 Merit Medical Europe • Horsterweg 24 • 6199 AC MAASTRICHT- AIRPORT • The Netherlands • +31 43 35 88 222
 Merit Medical Ireland Ltd. • Parkmore Business Park West • Galway, Ireland • +353 (0) 91 703 733
 Free phone from specific country:
 Belgium 0800 72906 • France 0800 91 60 30 • Germany 0800 182 0871 • Luxembourg 0800 22522 • Netherlands 0800 022 8184 • UK 0800 973 115



NEW!

8.3F KIT COMPONENTS:

- 8.3F (40cm) pigtail or straight catheter (with marker band)
- 8.5F (22cm) dilator
- .035" (0.89mm) x 80cm, double ended (J & straight) standard PTFE guide wire
- 25G x 1 1/2" (3.8cm) anesthesia needle
- 18G x 1 1/2" (3.8cm) medication needle
- 18G x 3.5" (9cm) access needle
- 18G x 5.9" (15cm) access needle
- #11 Futura safety scalpel
- 10ml luer lock syringe
- 20ml luer lock syringe
- 60ml luer lock syringe
- 1000ml disposal bag with drainage spout
- (2) Large bore three-way stopcocks
- Double-ended alligator clip cable
- (10) 4" x 4" (10cm x 10cm) gauze sponges
- 60cm x 60cm fenestrated drape
- 46cm x 60cm absorbent barrier tri-layer drape
- Suture 3-0 silk with curved cutting needle
- PAL label set – (4) yellow saline, (4) red specimen, (4) blue lidocaine
- (2) ChloroPrep® 3ml (orange tint) for use in PC801 and PC802



6F KIT COMPONENTS:

- 6F (60cm) high-flow pigtail or straight catheter
- 6F (22cm) dilator
- .035" (0.89mm)x 80cm, double ended (J and straight) heavy duty PTFE guidewire
- 22G x 2 3/8" (6cm) anesthesia needle
- 19G x 1 1/2" (3.8cm) medication needle
- 18G x 3.5" (9cm) access needle
- 18G x 5.9" (15cm) access needle
- #11 scalpel
- 10ml luer lock syringe
- 60ml luer lock syringe
- 1400ml disposal bag with connecting tube
- Three-way stopcock
- Double-ended alligator clip cable
- Ten 4" x 4" (10cm x 10cm) gauze sponges
- 40cm x 40cm fenestrated drape
- 46cm x 60cm absorbent barrier tri-layer drape

www.merit.com