

# Medela DCS



Medela's Disposable Collection System supports the safe collection and disposal of suctioned fluids.

## Disposable Collection System



Ref. 077.0082 (1.5l)  
Ref. 077.0085 (2.5l)

### Suction jar PC for disposable liners

- Reusable jars for holding disposable liners
- Made of unbreakable and easy to clean polycarbonate
- Sizes: 1.5 and 2.5 liter



Ref. 077.0083 (1.5l)  
Ref. 077.0086 (2.5l)

### Disposable suction liners

- Disposable solution for collecting secretions
- 2 layered foil with integrated overflow protection and lid
- Sizes: 1.5 and 2.5 liter



Ref. 077.0084 (1.5l)  
Ref. 077.0087 (2.5l)

### Disposable suction liners with solidifier

- Disposable solution for collecting secretions
- 2 layered foil with integrated overflow protection and lid
- Contains solidifier for leak-proof transport and disposal
- Sizes: 1.5 and 2.5 liter



Ref. 077.0094

### Disposable specimen cup

- For collecting tissue samples
- Can be attached to patient connection of lid
- Made of polypropylene
- Package of 10 pcs



Ref. 077.0188 A  
Ref. 077.0190 B

### Trolley

- For use with Medela Fluid Collection Systems
- Holds up to 5 suction jars/liners
- Mobile with 5 castors (with brakes)
- Dimensions: 85 x 60 x 60 cm