

H&H PriMed Gauze

Sterile



Item No: CL-0017
Manufacturer: H&H
Certificates: CE

Size: $7,62 \times 5,08 \times 2,54$ cm

Color: White Country of origin: China

NSN no: 6510 01 503 2117

Tariff: 40149050

Weight: 39g

H&H Primed Gauze takes care of moderate to severe wounds involving heavy blood loss.

Made with high quality cotton for maximum absorption and stability, H&H Primed Gauze is vacuum packed small for easy storage yet expands to cover large wound areas.

It is sterile and a perfect companion to our compression bandages to stop severe blood loss and save lives.

45 96 90 70 90 ap-services.dk Aulumvej 7 7550 Sørvad Denmarl