

TVAC SUCTION



Item No:	AW-000111
Manufacturer:	Fenton Pharmaceuticals
Certificates:	CE
Size:	H25cm x W28cm
Color:	N/A
Country of origin:	China
NSN no:	6515-01-541-7206
Tariff:	90189084
Weight:	190g

TVAC disposable portable suction unit with its resilient silicone bulb performs reliably in even the harshest climates without loss of suction.

Strong silicone bulb creates vacuum pressure of approx. 110mmHg. Collection bag volume >1,000ml.

Contoured suction tube tip prevents tissue damage during applications.