

PRODUCT DATA SHEET

INTROCAN SAFETY-W PUR 24G, 0.7X19MM-EU(4253523-01)



Attribute Name	Attribute Value
Brand name	Introcan Safety®
Product name (HC)	Introcan Safety®-W
Product type	IV catheter
Wing	Wing
Injection port	No
Flow rate	22.0 ml/min
Connection	Luer-Lock/female
Inner-ø catheter (metr.)	0.50 mm
Residual volume	0.093 ml
Color	yellow
Material of catheter	PUR
Safety clip	Yes
X-ray contrastable marking	X-ray contrast thread
Type of packaging	Blister
Product (unpack.) contains PVC	No
Product (unpack.) cont. Latex	No
Product (unpack.) cont. DEHP	No

Attribute Name	Attribute Value
Packaging contains PVC	No
Packaging contains DEHP	No
Packaging contains Latex	No
Product contains PVC	No
Product contains DEHP	No
Product contains Latex	No
Product contains Nickel	Yes
Sterility	sterile
Type of Sterilization	EO
UMDNS	10582-Cannulae, Venous
UMDNS-Version	1.0
Disposable / Reusable	disposable
Application	Vascular access, Peripheral venous access
eCl@ss Code (Rel. 11)	34220309
eCl@ss Description (Rel. 11)	Safety indw. venous cann.
Catheter-Gauge	24
Outer-ø catheter (metr.)	0.70 mm
Length of catheter	19.0 mm
Length of catheter (inch)	3/4"